

# PHILIPS

## sonicare



Sonic electric  
toothbrush with app

DiamondClean  
Smart 9300

Up to 20x more plaque removal\*

Up to 15x healthier gums\*\*

Up to 100% less stains\*\*\*

Visible pressure alert



HX9903/05

## Complete care for a healthier mouth

Feel confident with 100% coverage

The DiamondClean Smart provides a deep clean for a whiter, healthier smile. High-performance brush heads let you focus on all areas of your oral health and our Smart Sensor technology gives you personalized feedback and coaching.

### Personalized in-app coaching to perfect brushing

- Smart guidance for gentle care
- Helps you improve your brushing technique
- Guides you and provides instant feedback as you brush
- Takes you back to any missed spots for 100% coverage
- Start taking better care of your trouble areas

### Optimize your brushing

- Automatically selects the best mode for you
- Make sure you're getting the most from your brush head
- Choose from 4 modes and 3 intensity settings
- Premium charging glass and travel case

### Premium performance

- Up to 4x more surface contact\*\*\*\* for an effortless clean
- All-in-One brush head for complete care
- Sonic brushing: the masterclass in effective cleaning

# DiamondClean Smart 9300

Sonic electric toothbrush with app

HX9903/05


## Highlights

### 4 modes, 3 intensity settings




With DiamondClean Smart you get the full works. Our 3 intensity settings let you boost your clean, and our 4 modes cover all your brushing needs; Clean mode for exceptional daily cleaning, White+ to remove stains, Deep Clean+ for an invigorating deep clean, and Gum Health to gently clean along your gum line.

### TouchUp feature



If you happen to miss a spot while brushing your teeth, your app's TouchUp feature will show you. You can then go back for a second pass, and be confident you are getting a complete clean, every time.

### Focus areas and goal-setting



Any trouble areas your dentist has pointed out? We'll highlight them on your in-app 3D mouth map, and remind you to pay extra attention to these areas.

### Design with storage in mind



Charging a toothbrush has never been so effortless. Simply place your toothbrush into the sleek glass holder after brushing, and it instantly starts charging. Not only does your glass charge fit stylishly into your bathroom, it also doubles as a rinsing glass. Enjoy two weeks of regular use from a single full charge. Our premium travel case lets you store your toothbrush hygienically when you're on the go.

# DiamondClean Smart 9300

Sonic electric toothbrush with app

HX9903/05

## Specifications

### Connectivity

**Bluetooth® wireless technology:** Connected brushing app

### Power

**Voltage:** 110-220 V

### Technical specifications

**Battery:** Rechargeable

**Operating time (full to empty):** 14 days\*\*\*\*\*

**Battery type:** Lithium ION

### Design and finishing

**Color:** White

### Service

**Warranty:** 2-year limited warranty

### Compatibility

**Android compatibility:** Android phones with Bluetooth 4.0 and above, Search Sonicare on Google Play Store for detailed OS compatibility.

**iOS compatibility:** iPhone with Bluetooth 4.0 and above, Search Sonicare on App Store for detailed iOS compatibility.

### Ease of use

**Display:** Illuminated display

**Battery indicator:** Illuminated icon indicates battery life

**Handle:** Slim ergonomic design

**Handle compatibility:** Easy click-on brush heads

**Brush head recognition:** Syncs with optimal mode

**Replacement reminder:** To always ensure best results, reminder icon lights up

### Items included

**Handle:** 1 DiamondClean Smart

**Travel case:** Premium travel case

**Brush heads:** 1 A3 Premium All-in-One, 1 C3 Premium Plaque Control

**Glass charger:** 1

### Cleaning performance

**Speed:** Up to 62000 brush movement/min

**Health benefits:** Up to 15x healthier gums\*\*

**Performance:** Removes up to 20x more plaque\*

**Timer:** BrushPacer and SmarTimer

**Whitening benefits:** Up to 100% less stains\*\*\*

**Pressure feedback:** Ring lights up in purple, Vibrates handle to alert user

**TouchUp:** Ensures 100% coverage

### Modes

**3 intensities:** High, Medium, Low

**Clean:** For exceptional everyday clean

**Deep Clean+:** For an invigorating deep clean

**Gum Health:** Special attention to molars

**White+:** To remove surface stains

### BrushSync Mode Pairing

**C3 Premium Plaque Control:** Pairs with Clean mode

**A3 Premium All-in-One:** Pairs with Clean mode

### Smart sensor technology

**Pressure sensor:** Alerts when brushing too hard

**Location sensor:** Tracks & improves coverage

**Scrubbing sensor:** Guides to reduce scrubbing

**3D mouth map:** Highlights trouble areas

**TouchUp:** To address missed spots

**Position sensor:** Guides to areas you missed



\* vs. a manual toothbrush.

\*\* when used with A3 brush head vs. a manual toothbrush.

\*\*\* vs. a manual toothbrush. In a lab study, actual results may vary.

\*\*\*\* than a DiamondClean brush head.

\*\*\*\*\* based on two periods of two-minute brushings per day on standard mode