

# PHILIPS

## sonicare



Sonic electric  
toothbrush with app

DiamondClean  
Smart 9750

Up to 20x more plaque removal\*

Up to 15x healthier gums\*\*

Up to 100% less stains\*\*\*

Visible pressure alert



HX9954/74

## Complete care for a healthier mouth

Feel confident with 100% coverage

The DiamondClean Smart provides a deep clean for a whiter, healthier smile. High-performance brush heads let you focus on all areas of your oral health and our Smart Sensor technology gives you personalized feedback and coaching.

### Personalized in-app coaching to perfect brushing

- Helps you improve your brushing technique
- Guides you and provides instant feedback as you brush
- Takes you back to any missed spots for 100% coverage
- Alerts you if you're pressing too hard
- Start taking better care of your trouble areas

### Optimize your brushing

- Make sure you're getting the most from your brush head
- Choose from 5 modes and 3 intensity settings
- Automatically selects the best mode for you
- Easy and stylish charging - home and away

### Premium performance

- Enjoy long-lasting fresh breath
- Up to 7x healthier gums in just 2 weeks\*\*\*\*\*
- Up to 4x more surface contact\*\*\*\* for an effortless clean


# DiamondClean Smart 9750

Sonic electric toothbrush with app

HX9954/74


## Highlights

### TouchUp feature




If you happen to miss a spot while brushing your teeth, your app's TouchUp feature will show you. You can then go back for a second pass, and be confident you are getting a complete clean, every time.

### Focus areas and goal-setting



Any trouble areas your dentist has pointed out? We'll highlight them on your in-app 3D mouth map, and remind you to pay extra attention to these areas.


### Too much pressure?



You might not notice when you brush too hard, but your DiamondClean Smart will. If you use too much pressure, the light ring on the end of your

handle will flash. This is a gentle reminder to ease off the pressure, and let your brush head do the work. 7 out of 10 people found this feature helped them become a better brusher.

### Design with storage in mind



You'll enjoy two weeks of regular use from a single full charge. However, if you run low while traveling, your deluxe USB travel case doubles as a charger. Place your toothbrush in the case and either connect it to your laptop or a wall socket. As for charging at home - simply place your brush in your sleek charger glass - it fits stylishly into your bathroom and also doubles as a rinsing glass. Keep your brush heads clean and ordered in your matching brush head holder.

# DiamondClean Smart 9750

Sonic electric toothbrush with app

HX9954/74

## Specifications

### Connectivity

**Bluetooth® wireless technology:** Connected brushing app

### Power

**Voltage:** 110-220 V

### Technical specifications

**Battery:** Rechargeable

**Operating time (full to empty):** 14 days\*\*\*\*\*

**Battery type:** Lithium ION

### Design and finishing

**Color:** Lunar Blue

### Service

**Warranty:** 2-year limited warranty

### Compatibility

**Android compatibility:** Android phones with Bluetooth 4.0 and above, Search Sonicare on Google Play Store for detailed OS compatibility.

**iOS compatibility:** iPhone with Bluetooth 4.0 and above, Search Sonicare on App Store for detailed iOS compatibility.

### Ease of use

**Display:** Illuminated display

**Battery indicator:** Illuminated icon indicates battery life

**Handle:** Slim ergonomic design

**Handle compatibility:** Easy click-on brush heads

**Brush head recognition:** Syncs with optimal mode

**Replacement reminder:** To always ensure best results, reminder icon lights up

### Items included

**Handle:** 1 DiamondClean Smart

**Travel case:** USB charging travel case

**Brush head holder:** 1

**Brush heads:** 1 C3 Premium Plaque Control, 1 G3 Premium Gum Care, 1 TongueCare+ tongue brush, 1 A3 Premium All-in-One

**Glass charger:** 1

### Cleaning performance

**Speed:** Up to 62000 brush movement/min

**Health benefits:** Up to 15x healthier gums\*\*

**Performance:** Removes up to 20x more plaque\*

**Timer:** BrushPacer and SmarTimer

**Whitening benefits:** Up to 100% less stains\*\*\*

**Pressure feedback:** Ring lights up in purple, Vibrates handle to alert user

**TouchUp:** Ensures 100% coverage

### Modes

**3 intensities:** High, Medium, Low

**Clean:** For exceptional everyday clean

**Deep Clean+:** For an invigorating deep clean

**Gum Health:** Special attention to molars

**TongueCare+:** For long lasting clean breath

**White+:** To remove surface stains

### BrushSync Mode Pairing

**C3 Premium Plaque Control:** Pairs with Clean mode

**G3 Premium Gum Care brush head:** Pairs with Gum Health mode

**TongueCare+ tongue brush:** Pairs with TongueCare mode

**A3 Premium All-in-One:** Pairs with Clean mode

### Smart sensor technology

**Pressure sensor:** Alerts when brushing too hard

**Location sensor:** Tracks & improves coverage

**Scrubbing sensor:** Guides to reduce scrubbing

**3D mouth map:** Highlights trouble areas

**TouchUp:** To address missed spots



\* vs. a manual toothbrush.

\*\* when used with A3 brush head vs. a manual toothbrush.

\*\*\* vs. a manual toothbrush. In a lab study, actual results may vary.

\*\*\*\* than a DiamondClean brush head.

\*\*\*\*\* in Gum Health mode vs a manual toothbrush; measured by GBI.

\*\*\*\*\* based on two periods of two-minute brushings per day on standard mode.