

# PHILIPS

## sonicare



Power Toothbrush  
with SenseIQ

DiamondClean  
Prestige 9900

Our most advanced toothbrush

SenseIQ senses, adapts, cares

All-in-One brush head

AI-powered Sonicare app



HX9994

## Exceptional care personalized to you with SenseIQ

Philips Sonicare DiamondClean Prestige 9900 with SenseIQ technology is our most advanced power toothbrush. It enhances our proven sonic technology with intelligence, sensing and adapting for holistic care so you always get it right.

### Always get it right with SenseIQ

- Adapts to you in real time when you brush too hard
- Real-time feedback on your brushing
- SenseIQ technology senses, adapts and cares
- Senses your brushing style

### Personalized clean without compromise

- All-in-One brush head for complete care
- Personalized guidance from the Philips Sonicare app
- Sonic brushing: the masterclass in effective cleaning

### Premium accessories

- Easy charging with the minimalist stand and base
- Prestige travel case with USB for go-anywhere style

### Always get it right with Sense

### Designed around you

# DiamondClean Prestige 9900

Power Toothbrush with SenselQ

HX9994/11


## Highlights

### SenselQ technology



As you brush, SenselQ technology senses pressure, motion, coverage and more up to 100 times per second. It then adapts the intensity if you push too hard. The intelligence of SenselQ takes the guesswork out of brushing modes, so you can experience effortless care and better brushing over time.

### Premium All-in-One brush head



Clean without compromise with the A3 Premium All-in-One brush head. Our all-time best has angled bristles to help remove up to 20x more plaque even in hard-to-reach areas\*. Triangular tips work to remove up to 100% more stains in less than two days\*\*\*. And longer bristles clean deep for gums that are up to 15x healthier in just two weeks\*\*. All that without ever switching brush heads. Did you know: brush heads become less effective after 3 months of use? Our BrushSync™ feature can remind you when it's almost time for a replacement.

### Personalized guidance



The Philips Sonicare app is powered by artificial intelligence and works in perfect harmony with the toothbrush. Use it to brush with real-time guidance about pressure, motion, coverage, duration and frequency. View progress reports by day, week, month and year. Get personalized recommendations and actionable steps for better brushing. Meanwhile, auto-syncing keeps brushing data up-to-date in the app, even when you don't brush with the app by your side.

### Senses your brushing style



While you clean, Philips Sonicare 9900 Prestige uses sensors to detect the pressure you apply, cleaning motion and coverage – up to 100 times per second. The sensors also keep an eye on how long and how often you brush.

### Adapts to you in real time



Most of us apply too much pressure during brushing. That's why our toothbrush automatically adjusts intensity when you push too hard. And with the free Philips Sonicare app, you can experience tailormade guidance and personalized recommendations to take your brushing to the next level.

### Proven sonic brushing



Brush away plaque with sonic technology – our benchmark in effective and gentle cleaning. While you brush, up to 62,000 bristle movements and dynamic fluid action help clean even hard-to-reach areas, leaving your mouth feeling exceptionally fresh and clean.


# DiamondClean Prestige 9900

Power Toothbrush with SenselQ

HX9994/11

## Highlights

### Real-time feedback



You might not notice when you brush too hard or start scrubbing, but your toothbrush will. The light ring on the end of your handle will gently remind you to ease off the pressure or scrubbing motion.

### Prestige travel case with USB



The elegant and compact travel case can be slipped into any bag for carefree convenience while away from home. The durable case is wrapped in soft-touch vegan-friendly leather. And for extra ease, it includes an integrated USB-

C port, and a USB charging cable, so you can easily recharge the toothbrush without removing it from the case.

### Easy charging

When it is time to charge the toothbrush, simply place it on the charging base. Perfect for admiring your toothbrush's crisp, modern design.

# DiamondClean Prestige 9900

Power Toothbrush with SenseIQ

HX9994/11

## Specifications

### Connectivity

**Bluetooth® wireless technology:** Connected brushing app

### Power

**Power supply:** 110-220V

### Technical specifications

**Operating time (full to empty):** 2+ weeks\*\*\*\*

**Battery:** Rechargeable

**Battery type:** Lithium ION

### Design and finishing

**Color:** Pink

### Service

**Warranty:** 2-year limited warranty

### Compatibility

**Android compatibility:** Android OS 8.0 or above

**iOS compatibility:** iPhones on iOS 13.0 or above

### Ease of use

**Handle:** Sleek and compact design, Seamless touch control

**Handle compatibility:** Easy click-on brush heads

**Battery indicator:** Illuminated icon indicator

**Travel case:** Elegant, compact charging case, Vegan-friendly leather

**BrushSync Replacement Reminder:** Reminder icon lights up, to always ensure best results

### Items included

**Handle:** 1 DiamondClean Prestige 9900

**Brush heads:** 1 A3 Premium All-in-One

**Prestige USB charging case:** 1

**Charger base:** 1

**USB cable:** 1

### Cleaning performance

**Health benefits:** Up to 15x healthier gums\*\*

**Performance:** Removes up to 20x more plaque\*

**Speed:** Up to 62000 brush movement/min

**Whitening benefits:** Up to 100% less stains\*\*\*

**Timer:** BrushPacer and SmarTimer

**Pressure feedback:** Ring lights up in purple, Vibrates handle to alert user

**Scrubbing feedback:** Ring lights up in amber, Vibrates handle to alert user, Option to turn on/off in app

### SenseIQ technology

**Pressure sensor:** Alerts when brushing too hard

**Position sensor:** Tracks brushing coverage, for each brushing session

**Scrubbing sensor:** Guides to reduce scrubbing

**3D mouth map:** Highlights missing areas

**Adaptive intensity:** Adjusts intensity level, when too much pressure applied

**Artificial Intelligence:** Detects and syncs brushing

### Sonicare App with A.I.

**Real time guidance:** Personalized recommendations, Day/month/year progress report

**Customizable modes:** Clean, White+, Gum health, Deep Clean, Sensitive

