

# PHILIPS

## sonicare



Patient

Demonstration Brush

In-Mouth Trial Unit

Up to 20 x more plaque removal

2 modes, 3 intensities

Visible pressure alert

21-day battery life



HX7408/01

## Gentle clean, even for sensitive teeth and gums

Up to 20 x more plaque removal\*

Brush for demonstrating the power of Sonicare, while in-surgery. Contains: 1 x Sonicare Series 6100 30 x A3 Brush Heads 30 x S2 Brush Heads

### Gentle, superior cleaning

- Complete care for a healthier mouth

### Brushing made effortless

- Always know when it's time for a new brush head
- Personalise your clean with 6 brushing settings
- Protect gums with the visible pressure sensor

### Superior clean

- Next-Generation Sonicare Technology

### Designed for your lifestyle

- Laser-welded for strength and durability
- Up to 3 weeks of brushing with one full charge
- Optimise brushing with Smartimer and BrushPacer

# In-Mouth Trial Unit

Patient Demonstration Brush

HX7408/01

## Highlights

### A3 All-in-One and S2 Sensitive



The A3 Premium All-in-One brush head is designed with patients in mind. With the A3 brush head, patients can achieve up to 20 x more plaque removal\*, up to 15 x healthier gums\*\* and remove up to 100% more stains\*\*\*. Also included: Thanks to the new S2 Sensitive Brush Head made using Sonicare's softest bristles, patients can achieve a gentle clean that removes up to 10 x more plaque in hard-to-reach areas\*.

### Next-Gen Sonicare technology



Enjoy a delightful brushing experience and new level of oral care with Next-Generation Sonicare technology. While you brush, the toothbrush's innovative magnetic system adjusts the brushing power during each cleaning session, delivering optimal performance of 62,000 bristle movements per minute, even in hard-to-reach areas. Meanwhile, dynamic fluid action supports the bristles to clean plaque, helping you achieve a complete clean every time you brush.


### 6 brushing settings



Designed for sensitive teeth and gums, the Philips Sonicare 6100 allows patients to choose between two brushing modes. Sensitive mode

provides gentle, thorough cleaning, while Clean mode focuses on plaque removal for a superior clean. For a personalised clean, patients can combine any mode with one of three intensity levels, from High to Low, for a brushing session that's just right.

### Light-ring pressure sensor




Patients may not notice when they're brushing too hard, but this toothbrush will. The smart optic sensor precisely detects excessive brushing pressure, then flashes purple while reducing vibration. Patients will always know when to adjust their brushing, and their gums will stay protected.

### BrushSync replacement reminder




Replacing your brush head regularly is an important part of maintaining good oral health. With BrushSync technology, it's easy for patients to keep their brush head performing at its best. By tracking usage frequency and brushing intensity, the innovative technology then reminds patients when it's time for a new brush head.

### 2-minute timer and brushing pacer




Two minutes is all patients need for a complete clean, so our Smartimer lets patients know when time's up. For an even clean throughout the mouth, the BrushPacer divides brushing sessions into segments, with a slight pause in vibration being a sign to move on. This supports patients in achieving an even clean every time they brush.

### Designed to last



We've re-designed our Philips Sonicare toothbrushes from the inside out to be strong, lightweight and long-lasting. Thanks to our advanced laser-welding process, our toothbrushes are designed to brush patients' teeth for years to come.

### Up to 3 weeks between charges



This toothbrush uses a high-performance battery that is engineered to enhance your brushing experience. With a single charge, patients can enjoy up to three weeks of regular brushing before needing to recharge it.

# In-Mouth Trial Unit

Patient Demonstration Brush

HX7408/01

## Specifications

### Intensities

**High:** To boost your clean

**Medium:** For an everyday clean

**Low:** For sensitive teeth and gums

### Smart sensor technology

**Replacement Reminder:** BrushSync lets you know when to replace your brush head

### Smart sensor technology

**Pressure feedback:** Light ring illuminates purple, Reduced vibration and pulsing sound

### Power

**Voltage:** 100-240 V

### Technical specifications

**Operating time (full to empty):** 21 days

**Battery:** Rechargeable

**Battery type:** Lithium ION

### Design and finishing

**Colour:** White

### Service

**Warranty:** 2-year limited warranty

### Ease of use

**Handle compatibility:** Easy click-on brush heads

**Handle:** Slim ergonomic design

**Battery indicator:** Illuminated icon indicates battery life

**Timer:** SmarTimer and BrushPacer for complete and thorough brushing

**USB-A Charging:** Power adapter not included

### Items included

**Charger:** 1

**Handle:** 1 6100 rechargeable toothbrush

**Brush head:** 30 S2 Sensitive, 30 A3 Premium All-in-One

### Cleaning performance

**Plaque removal:** Up to 20 x more effective\*

**Speed:** 62,000 brush movements/min

**Gum Health:** Up to 15 x healthier gums\*\*\*

**Whitening:** Up to 100% fewer stains\*\*

### Modes

**Clean:** For exceptional everyday cleaning

**Sensitive:** For sensitive teeth and gums

