

PHILIPS

Refurbished IPL
Hair removal device

Lumea IPL 8000
Series

SmartSkin sensor

4 intelligent attachments: body,
face, bikini, underarms

Lumea IPL App

Corded use



BRI948/00R1



Tailored to you with SenseIQ

Enjoy 18 months of visibly hair-free smooth skin*

Enjoy a quick and personalised treatment with our Lumea IPL 8000 Series. SenseIQ technology with smart attachments and Lumea IPL app for long-lasting smooth skin.

Gentle and effective for long-lasting smooth skin

- Treat only 2 x a month for fast results
- Developed with dermatologists to be easy and effective
- Extra-long cable for added flexibility while treating

Full solution for face and body

- Full solution for face and body with 4 smart attachments
- Suitable for most skin tones and hair colours

Refurbished Lumea IPL

- Same results, better for the environment

Personalised IPL hair removal treatment

- Gentle and comfortable with SenseIQ
- Optimise your routine with the Philips Lumea IPL app


Highlights

Refurbished Lumea IPL



We love to reuse returned products and give them a second life, because the materials and parts we use to make them have the biggest CO2 impact. Choosing refurbished instead of brand new reduces the carbon footprint by 78%.***

2 x a month for fast results



Every 2 weeks to start with — that's half as many treatments as other brands. Followed by touch-ups just once a month. That's it. Covering both lower legs takes 8.5 minutes.

SenseIQ technology



Lumea 8000 Series has five easily adjustable light settings. Our SmartSkin sensor reads your skin tone and helps you find the most comfortable setting. Intelligent attachments adapt treatment for each body area.

Extra-long cable



Convenient to use thanks to an extra long cable for easy access and enhanced maneuverability.

4 smart attachments



Specially designed attachments perfectly fit the curves of your body and trigger the most effective programmes for each body area when they are connected. Face: Flat design and small window with UV filter. Body: Curved inward with a large window. Bikini and Armpit: Curved outward for hard-to-reach areas.


Developed with dermatologists



As a leader in health technology, Philips developed Lumea IPL in consultation with dermatologists, for easy and effective use in the safety of your home. Derived from

technology used in professional salons, Lumea IPL offers gentle treatment, even on sensitive areas.

Philips Lumea IPL app



Our free coaching app helps you plan and stick to your treatment schedule, then takes you through each session step by step. Downloaded by more than 2.1 million users.

For most skin and hair types



IPL needs contrast between the pigment in the hair colour and the pigment in the skin tone, it therefore works on naturally dark blonde, brown and black hair and on skin tones from fair to dark brown (I-V).

Specifications

Technical specs. attachments

Body treatment attachment: Shape: Convex curved, Window size: 4.1 cm², Tailored treatment for body: legs, arms and stomach

Face treatment attachment: Shape: Flat, Window size: 2 cm², extra filter, Tailored treatment for face: upper lip, chin and jawline

Bikini treatment attachment: Shape: Concave curved, Window size: 3 cm², Tailored treatment for sensitive bikini area

Underarm treatment attachment: Shape: Concave curved, Window size: 3 cm², Tailored treatment for underarms

Safety and adjustable settings

Integrated UV filter: Protects skin from UV light

5 light energy settings: Adjustable to your individual skin tone

Integrated safety system: Prevents unintentional flashing

SmartSkin sensor: The right setting on demand

Application mode

Slide and Flash: For quick application

Stamp and Flash: For treatment on small areas

Corded/cordless use: Corded use

Technical specifications

Voltage: 100–240 V

High performance lamp: Built to last, 450,000 flashes, equivalent to 39 years of lamp lifetime**

Service

Warranty: 2 year global warranty + 3 extra years of warranty upon product registration within 90 days

Items included

Instructions for use: User manual

Storage: Beauty pouch

Adapter: 19.5 V / 4000 mA

Application time

Underarms: 2.5 min

Bikini line: 2 min

Face areas: 1.5 min

Lower legs: 8.5 min



** After 3 treatments. Median result 86% hair reduction on lower legs at 18 months. Individual results may vary. Based on a survey of 90 consumers.

** When following the treatment schedule. Calculated for use on lower legs, bikini, armpits and face. The lamp lifetime does not extend the Philips 2-year worldwide guarantee

*** In scope: product, packaging materials and Philips reverse logistics of Lumea S8000 returns. Out of scope: chemicals, energy needed for reprocessing, inbound part and material transport. Calculation based on ReCiPe2016 and ecoinvent3.6 database, assuming all of the accessories, the adapter and the packaging need replacements.