

PHILIPS

sonicare

Standard sonic
toothbrush heads

W2 Optimal White

5-pack

Standard size

Click-on

BrushSync mode pairing



HX6065/10

Up to 100% whiter teeth in just one week*

Advanced cleaning for stain removal and whiter teeth

The W2 Optimal White brush head is perfect for those who want to go beyond deep cleaning to remove surface stains for a white, radiant smile. This brush head is also great for maintaining brightness between professional whitening treatments

Innovative technology

- Automatically selects optimal mode for ultimate results**
- Philips Sonicare's advanced sonic technology

Designed to optimise performance

- Tested to meet your oral health needs
- Click-on design for simple brush head placement
- Always know when to replace. Always get an effective clean.

For exceptionally whiter teeth

- Diamond-shaped bristles give you 100% whiter teeth in 1 week

For exceptional care of teeth and gums

- Up to 7x more plaque removal than a manual toothbrush

W2 Optimal White

Standard sonic toothbrush heads

HX6065/10


Highlights

Say goodbye to plaque



Densely-packed, high-quality bristles remove up to 7x more plaque than a manual toothbrush.

Whiter teeth, fast



The DiamondClean's stain-removal pad, made from densely-packed diamond-shaped bristles, removes surface stains caused by food and drink. You'll notice a whiter smile, up to 100%*, in just 7 days.

BrushSync™ mode-pairing



You'll always get the best possible clean with our BrushSync™ mode-pairing feature**. The Philips Sonicare W2 Optimal White brush head syncs

with your BrushSync™ enabled Philips Sonicare toothbrush handle**, selecting the optimal brushing mode and intensity level for exceptional whitening. All you need to do is start brushing.

Brush head replacements



Brush heads become less effective after 3 months' use, but with BrushSync™ you'll be reminded before this happens. Your smart toothbrush will track how often and how hard you brush and will notify you when it's time for a replacement. Don't have a smart Philips Sonicare toothbrush? Just monitor your blue replacement bristles — when they turn white you'll know it's time for a fresh brush head.

Brushing with peace of mind



All Philips Sonicare brush heads are safe and gentle on teeth and gums. During tests, each brush head is put through its paces to ensure durability and an exceptional performance every time you brush.

Philips Sonicare technology



Philips Sonicare's advanced sonic technology pulses water between your teeth, and its brush strokes break plaque up and sweep it away for an exceptional daily clean.

Easy to click on



Your W DiamondClean brush head fits perfectly with any Philips Sonicare toothbrush handle, except Philips One and Essence. Simply click on and off for easy replacement and cleaning.

W2 Optimal White

Standard sonic toothbrush heads

HX6065/10

Specifications

Consumer Trade Item

Height: 4.4 cm

Width: 6 cm

Length: 22.1 cm

Net weight: 0.036 kg

Gross weight: 0.056 kg

GTIN: 8710103911906

Country of origin: DE

Harmonized system code: 96032100

Outer Carton

Height: 7.5 cm

Width: 22 cm

Length: 38 cm

Net weight: 0.36 kg

Gross weight: 0.689 kg

GTIN: 18710103911903

Design and finishing

Colour: White

Bristle stiffness feel: Medium

Reminder bristles: Blue bristle colour fades away

Size: Standard

Pallet Layer

Height: 7.5 cm

Width: 80 cm

Length: 120 cm

Compatibility

Brush head system: Click-on

Suitable for these models: 2 Series plaque defence, 2 Series plaque defence, 3 Series gum health, DiamondClean, DiamondClean Smart, EasyClean, Essence+, FlexCare, FlexCare Platinum,

FlexCare Platinum Connected, FlexCare+, for Kids, HealthyWhite, HealthyWhite+, PowerUp, ProtectiveClean

[BrushSync mode pairing](#)

Not suitable for: Philips One

Items included

Brush heads: 5 W2 Optimal White

Quality and performance

Replacement: Every 3 months

Tested: for optimal usage

Health benefits

[Gum health](#): Helps improve gum health

[Plaque removal](#): Removes up to 7x more plaque*

[Whitening](#): Whiter teeth in just one week

