**PHILIPS**

**sonicare**

— |

| —

— |

| —

English 6  
Azərbaycanca 9  
Čeština 12  
Dansk 15  
Deutsch 18  
Eesti 22  
Español 25  
Français 28  
Hrvatski 31  
Italiano 34  
Latviešu 37  
Lietuviškai 40  
Magyar 43  
Nederlands 46  
Norsk 49  
O'zbek 52  
Polski 55  
Português 59  
Română 62  
Shqip 65  
Slovenščina 68  
Slovensky 71  
Srpski 74  
Suomi 77  
Svenska 80  
Türkmençe 83  
Ελληνικά 86  
Български 89  
Кыргызча 92  
Русский 95  
Тоҷикӣ 99  
Українська 102  
Қазақша 105

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1

## Important safety information

Only use the product for its intended purpose. Read this important information carefully before you use the product and its batteries and accessories, and save it for future reference. Misuse can lead to hazards or serious injuries.

Note: The accessories and features of different models vary. For example, some models come with a USB wall adapter or are equipped with multiple modes.

## Warnings

- Keep the charger away from water.
- After cleaning, make sure the charger is completely dry before you connect it to the wall socket.
- This appliance can be used by children and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved. Cleaning and user maintenance shall not be made by children without supervision.
- Children shall not play with the appliance.
- Only use original Philips accessories or consumables. Only use the USB charger provided with the product.
- Do not charge the product outdoors or near heated surfaces.
- Do not clean any part of the product in the dishwasher.
- To avoid physical injury, do not fill the reservoir with hot water. Do not put foreign objects in the reservoir.
- The maximum altitude of use is 4420 meters.
- This product has only been designed for cleaning teeth and gums. Do not use it for any other purpose. (for e.g. eyes, nose and ears).
- Stop using the appliance and consult your dentist/doctor if excessive bleeding occurs after use, if bleeding continues to occur after 1 week of use or if you experience discomfort or pain.
- If you have had oral or gum surgery in the previous 2 months, consult your dentist before you use this appliance.
- If you have a pacemaker or other implanted device, contact your doctor or the manufacturer of the implanted device prior to use.
- If you have medical concerns, consult your doctor before you use this appliance.
- This appliance is a personal care device and is not intended for use on multiple patients in a dental practice or institution.
- Charge and use the product at a temperature between 5 °C and 40 °C (41 °F and 104 °F).

## Battery safety instructions

- This appliance contains batteries that are only to be removed by skilled persons
- Only use the product for its intended purpose. Read this information carefully before you use the product and its batteries and accessories, and save it for future reference. Misuse can lead to hazards or serious injuries. Accessories supplied may vary for different products.
- Only use original Philips accessories and consumables. Only use detachable supply unit WAAxxxx.
- Keep product and batteries away from fire and do not expose them to direct sunlight or high temperatures.
- If the product becomes abnormally hot or smelly, changes color or if charging takes longer than usual, stop using and charging the product and contact Philips.
- Do not place products and their batteries in microwave ovens or on induction cookers.
- Do not open, modify, pierce, damage or dismantle the product or battery to prevent batteries from heating up or releasing toxic or hazardous substances. Do not short-circuit, overcharge or reverse charge batteries.
- This appliance contains batteries that are non-replaceable. When the battery is at end of life, the appliance shall be properly disposed of, see section on recycling.
- To avoid accidental short-circuiting of batteries after removal, do not let battery terminals come into contact with metal object (e.g. coins, hairpins, rings). Do not wrap batteries in aluminium foil. Tape battery terminals or put batteries in a plastic bag before you discard them.

## Electromagnetic fields (EMF)

This Philips appliance complies with all applicable standards and regulations regarding exposure to electromagnetic fields.

## Warranty and support

For product information, support and online user instructions, please visit: [www.philips.com/support](http://www.philips.com/support) or read the international warranty leaflet.

## Warranty restrictions

The terms of the international warranty do not cover the following:

- Reservoir
- Nozzles
- Damage caused by use of unauthorized replacement parts or unauthorized nozzles

## 8 English

- Damage caused by misuse, abuse, neglect, alterations or unauthorized repair
- Normal wear and tear, including chips, scratches, abrasions, discoloration or fading

### Recycling

- This symbol means that electrical products and batteries shall not be disposed of with normal household waste (Fig. 1).
- Follow your country's rules for the separate collection of electrical products and batteries.

CONTAINS LITHIUM-ION RECHARGEABLE CELL. MUST BE RECYCLED OR DISPOSED OF PROPERLY.

The built-in rechargeable battery must be removed by a qualified professional when the product is discarded. Instructions for removal of built-in rechargeable batteries can be found on [www.philips.com/support](http://www.philips.com/support).

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©2022 Koninklijke Philips N.V. (KPNV). All rights reserved. Philips and the Philips shield are trademarks of KPNV. Sonicare and the Sonicare logo are trademarks of Philips Oral Healthcare, LLC and/or KPNV.