

# PHILIPS

## sonicare

Sonic electric  
toothbrush with app

DiamondClean  
Smart 9400

Personalised coaching

Smart brush head sensors

Smart brush head recognition

4 modes, 3 intensities

HX9917/90



## Complete care for a healthier mouth

Feel confident, with 100% coverage

The DiamondClean Smart 9400 is our best toothbrush, for complete care. High-performance brush heads let you focus on all areas of your oral health and our Smart Sensor technology gives you personalised feedback and coaching.

### Proven to improve oral health

- Removes up to 10 x more plaque\* for a deep clean

### Personalised in-app coaching to perfect brushing

- Guides your brushing and highlights missed areas
- Guides you to reduce your scrubbing
- Takes you back to any missed spots for 100% coverage
- Start taking better care of your trouble areas
- Alerts you if you're pressing too hard

### Optimise your brushing

- Make sure you're getting the most from your brush head
- Easy and stylish charging — home and away
- Brush heads automatically select the optimal settings
- Choose from 4 modes and 3 intensity settings

# DiamondClean Smart 9400

Sonic electric toothbrush with app

HX9917/90

## Highlights

### Say goodbye to plaque




Click on the G3 Premium Plaque Control brush head to experience our deepest clean. With soft flexible sides, the bristles contour to the shape of each tooth's surface to provide 4 x more surface contact\*\* and remove up to 10 x more plaque from hard-to-reach areas\*.

### 4 modes, 3 intensity settings




With DiamondClean Smart you get the full works. Our 3 intensity settings let you boost your clean and our 4 modes cover all your brushing needs; Clean mode for exceptional daily cleaning, White+ to remove stains, Deep Clean+ for an invigorating deep clean, and Gum Health to gently clean along your gum line.

### Never miss a spot




The Philips Sonicare app lets you know when you've achieved a thorough clean, and coaches you to be a better, more mindful brusher.

### TouchUp feature



If you happen to miss a spot while brushing your teeth, your app's TouchUp feature will show you. You can then go back for a second pass, and be confident you are getting a complete clean, every time.


### Smart replacement reminders



Brush heads wear out over time and become less effective at removing plaque. With smart brush head recognition you'll never lose track of when


to replace your brush head. Your toothbrush tracks how long and hard you've brushed with each brush head, and will alert you when it's time to replace. You can also monitor the lifetime of your brush heads in the Philips Sonicare app, and even order replacements so you're never caught short.

### Focus areas and goal-setting



Any trouble areas your dentist has pointed out? We'll highlight them on your in-app 3D mouth map and remind you to pay extra attention to these areas.

### Too much pressure?



You might not notice when you brush too hard, but your DiamondClean Smart will. If you use too much pressure, the light ring on the end of your handle will flash. This is a gentle reminder to ease off the pressure and let your brush head do the work. 7 out of 10 people found that this feature helped them to become a better brusher.

# DiamondClean Smart 9400

Sonic electric toothbrush with app


HX9917/90

## Highlights

### Easy does it

With an electric toothbrush, you let the brush do the work. We've built a scrubbing sensor into your handle as a gentle reminder to reduce scrubbing. This way, you can improve your technique and get a gentler, more effective clean.


### BrushSync mode pairing



Wondering which mode and intensity you should use? Wonder no longer. A microchip in your brush head tells your DiamondClean Smart

which one you're using. So if you click on a gum-care brush head, for example, your toothbrush knows to select the optimal mode and intensity to provide gentle yet effective care for your gums. All you need to do is to press the power button.

### Designed with storage in mind



You'll enjoy two weeks of regular use from a single full charge. However, if you run low while travelling, your deluxe USB travel case doubles as

a charger. Place your toothbrush in the case and either connect it to your laptop or a wall socket. As for charging at home, simply place your brush in your sleek charger glass—it fits stylishly into your bathroom and also doubles as a rinsing glass.

# DiamondClean Smart 9400

Sonic electric toothbrush with app

HX9917/90

## Specifications

### Connectivity

**Bluetooth® wireless technology:** Connected brushing app

### Power

**Voltage:** 110-220 V

### Technical specifications

**Battery:** Rechargeable

**Operating time (full to empty):** 14 days\*\*\*

**Energy consumption:** Standby without display  
<0.5 W

**Battery type:** Lithium ION

### Design and finishing

**Colour:** Grey

### Service

**Warranty:** 2-year limited warranty

### Compatibility

**Android compatibility:** Android OS 8.0 or above

**iOS compatibility:** iPhones on iOS 13.0 or above

### Ease of use

**Handle:** Slim ergonomic design

**Battery indicator:** Illuminated icon indicates battery life

**Handle compatibility:** Easy click-on brush heads

**Brush head recognition:** Syncs with optimal mode

**Replacement reminder:** To always ensure best results, reminder icon lights up

### Items included

**Handle:** 1 DiamondClean Smart

**Brush heads:** 1 All-in-1 Premium Brush Head, 1 G3 Premium Gum Care

**Travel case:** USB charging travel case

**Glass charger:** 1

### Cleaning performance

**Speed:** Up to 62,000 brush movements/min

**Performance:** Removes up to 10 x more plaque\*

**Timer:** BrushPacer and SmarTimer

**Pressure feedback:** Ring lights up in purple, Vibrates handle to alert user

**TouchUp:** Ensures 100% coverage

### Modes

**Clean:** For exceptional everyday cleaning

**3 intensities:** Low, Medium, High

**Deep Clean+:** For an invigorating deep clean

**Gum Health:** For a gentle yet effective clean for gums

**White+:** To remove surface stains

### BrushSync Mode Pairing

**G3 Premium Gum Care brush head:** Pairs with Gum Health mode

**All-in-One Premium:** Pairs with clean mode

### Smart sensor technology

**Pressure sensor:** Alerts when brushing too hard

**Location sensor:** Tracks and improves coverage

**Scrubbing sensor:** Guides to reduce scrubbing

**3D mouth map:** Highlights trouble areas

**TouchUp:** To address missed spots

### Software Support

**Software updates:** Philips offers relevant software updates for a period of 2 years after the date of purchase.

