

PHILIPS

sonicare



Power Toothbrush
with SenseIQ

DiamondClean
Prestige 9900

Our most advanced toothbrush
SenseIQ senses, adapts, cares
All-in-One brush head
AI-powered Sonicare app



HX9992/45

Exceptional care personalised to you with SenseIQ

Philips Sonicare DiamondClean Prestige 9900 with SenseIQ technology is our most advanced power toothbrush. It enhances our proven sonic technology with intelligence, sensing and adapting for holistic care so you always get it right.

Designed around you

- Sleek and compact design

Always get it right with SenseIQ

- Adapts to you in real time when you brush too hard
- Cares intuitively for teeth and gums
- Real-time feedback on your brushing
- Senses your brushing style
- SenseIQ technology senses, adapts and cares

Personalised clean without compromise

- All-in-One brush head for complete care
- Personalised guidance from the Philips Sonicare app
- Sonic brushing: the masterclass in effective cleaning

Premium accessories

- Easy charging with the minimalist stand and base
- Prestige travel case with USB for go-anywhere style


DiamondClean Prestige 9900

Power Toothbrush with SenseIQ

HX9992/45


Highlights

SenseIQ technology



As you brush, SenseIQ technology senses pressure, motion, coverage and more up to 100 times per second. It then adapts the intensity if you push too hard. The intelligence of SenseIQ takes the guesswork out of brushing modes, so you can experience effortless care and better brushing over time.

Premium All-in-One brush head



Clean without compromise with the A3 Premium All-in-One toothbrush head. Our all-time best has angled bristles to help remove up to 20 x more plaque even in hard-to-reach areas*. Triangular tips work to remove up to 100% more stains in less than two days***. And longer bristles clean deep for gums that are up to 15 x healthier in just two weeks**. All that without ever switching toothbrush heads. Did you know: toothbrush heads become less effective after 3 months of use? Our BrushSync™ feature can remind you when it's almost time for a replacement.

Personalised guidance



The Philips Sonicare app is powered by artificial intelligence and works in perfect harmony with the toothbrush. Use it to brush with real-time guidance about pressure, motion, coverage, duration and frequency. View progress reports by day, week, month and year. Get personalised recommendations and actionable steps for better brushing. Meanwhile, auto-syncing keeps brushing data up to date in the app, even when you don't brush with the app by your side.

Senses your brushing style



While you clean, Philips Sonicare DiamondClean Prestige 9900 uses sensors to detect the pressure you apply, cleaning motion and coverage up to 100 times per second. The sensors also keep an eye on how long and how often you brush.

Adapts to you in real time



Most of us apply too much pressure during brushing. That's why our toothbrush automatically adjusts intensity when you push too hard. And with the free Philips Sonicare app,

you can experience tailor-made guidance and personalised recommendations to take your brushing to the next level.

Cares for your oral health




Philips Sonicare DiamondClean Prestige 9900 electric toothbrush is designed around you. Its state-of-the-art technology is so in tune with you and effortless to use that you barely know it is there. Luckily, the results say it all.

Proven sonic brushing



Brush away plaque with sonic technology – our benchmark in effective and gentle cleaning. While you brush, up to 62,000 bristle movements and dynamic fluid action help clean even hard-to-reach areas, leaving your mouth feeling exceptionally fresh and clean.

Real-time feedback



You might not notice when you brush too hard or start scrubbing, but your toothbrush will. The light ring on the end of your handle will gently remind you to ease off the pressure or scrubbing motion.

DiamondClean Prestige 9900

Power Toothbrush with SenseIQ

HX9992/45

Highlights

Prestige travel case with USB



The elegant and compact travel case can be slipped into any bag for care-free convenience while away from home. The durable case is wrapped in soft-touch, vegan-friendly leather. And for extra ease, it includes an integrated USB-

C port and a USB charging cable, so you can easily recharge the toothbrush without removing it from the case.

Easy charging



When it is time to charge the toothbrush, simply place it on the charging base. Perfect for admiring your toothbrush's crisp, modern design.

Sleek and compact



The Sonicare 9900 Prestige experience is about big technology delivered with style and compactness. The toothbrush is 18 mm shorter than DiamondClean and features premium materials to make this luxurious device easy on the eye and pleasant to the touch.

DiamondClean Prestige 9900

Power Toothbrush with SenseIQ

HX9992/45

Specifications

Connectivity

Bluetooth® wireless technology: Connected brushing app

Power

Power supply: 110-220 V

Technical specifications

Operating time (full to empty): 2+ weeks****

Battery: Rechargeable

Battery type: Lithium ION

Energy consumption: Standby without display
<0.5 W

Design and finishing

Colour: Midnight Blue

Service

Warranty: 2-year limited warranty

Compatibility

Android compatibility: Android OS 8.0 or above

iOS compatibility: iPhones on iOS 13.0 or above

Ease of use

Handle compatibility: Easy click-on brush heads

Handle: Sleek and compact design, Seamless touch control

Battery indicator: Illuminated icon indicator

Travel case: Elegant, compact charging case, Vegan-friendly leather

BrushSync Replacement Reminder: Reminder icon lights up, to always ensure best results

Items included

Brush heads: 4 A3 Premium All-in-One

Handle: 1 DiamondClean Prestige 9900

Prestige USB charging case: 1

Charger base: 1

USB cable: 1

Cleaning performance

Performance: Removes up to 20 x more plaque*

Speed: Up to 62,000 brush movements/min

Health benefits: Up to 15 x healthier gums**

Timer: BrushPacer and SmarTimer

Whitening benefits: Up to 100% less stains***

Pressure feedback: Ring lights up in purple, Vibrates handle to alert user

Scrubbing feedback: Ring lights up in amber, Vibrates handle to alert user, Option to turn on/off in app

Software Support

Software updates: Philips offers relevant software updates for a period of 2 years after the date of purchase.

SenseIQ technology

Pressure sensor: Alerts when brushing too hard

Position sensor: Tracks brushing coverage, for each brushing session

Scrubbing sensor: Guides to reduce scrubbing

3D mouth map: Highlights missing areas

Adaptive intensity: Adjusts intensity level, when too much pressure applied

Artificial Intelligence: Detects and syncs brushing

Sonicare App with A.I.

Real-time guidance: Personalised recommendations, Day/month/year progress report

Customisable modes: Clean, White+, Gum health, Deep Clean, Sensitive

