



PHILIPS
sonicare

Healthier teeth &
gums **in 2 weeks**

Powerfully effective, yet remarkably gentle



Performs up to
62,000 brush
movements
per minute



Removes
up to 10x
more plaque
between teeth



The #1 most-
recommended
sonic power
toothbrush
by dental
professionals
worldwide



Noticeably
whiter teeth
from day 1





Feel the sonic difference

Innovative sonic technology pulses water between your teeth

Philips Sonicare's advanced sonic technology pulses water between teeth, and its brush strokes break up plaque and sweep it away for that 'just been to the hygienist' professional feeling, every time.

Sonicare solutions for you at any price

There is powerful sonic technology in all of the handles in our range, plus any Sonicare brush head will fit on any Sonicare handle. We're proud that there is a Sonicare solution to suit any budget.



Sonicare For Kids Connected

The parent's choice
for better brushing

- 98% of parents saw improved brushing habits²
- Built-in Bluetooth and interactive coaching app to help kids learn about oral care through play
- Two minute timer
- Design your own brush with fun, changeable stickers

2. Philips Sonicare in a home survey of US dental professionals with children aged 4–10



Series 3100

Goodbye manual toothbrush.
Hello Sonic technology.

- Helps reduce cavities
- SmartTimer and QuadPacer
- Brush head replacement reminder
- Two week battery life
- Goodbye adapter.
Hello USB charging
- Built-in pressure sensor







Series 5300

Superior plaque removal
and whiter teeth

- 7x more plaque removal
- Adaptive brushing technology
- 2 intensities
- SmartTimer and SextantPacer
- Optical pressure sensor
- 21 day battery life
- Brush head replacement reminder



Series 6100

Gentle clean, even for
sensitive teeth and gums

- 10x more plaque removal
- Adaptive brushing technology
- 2 modes, 3 intensities
- Built-in pressure sensor with visual and haptic response
- SmartTimer and SextantPacer
- 21 day battery life
- Travel case



Series 7100

Perfect clean with up to 7x healthier gums in just 2 weeks

- 10x more plaque removal
- Adaptive brushing technology
- 4 modes and 3 intensities for the most personalised clean
- Built-in pressure sensor with visual and haptic response
- SmarTimer and SextantPacer
- Connected bluetooth brushing app
- 21 day battery life
- Charging travel case




Cordless Power Flosser 3000

Enjoy an effortless, thorough clean between your teeth and improve your gum health.

- Remove 99% of plaque in treated areas⁴
- Access hard-to-reach areas with a nozzle that rotates 360 degrees
- Quad Stream nozzle for maximum coverage
- Easy to fill 250ml reservoir for a full 60 second clean
- 2-week battery life

4. In an in-vitro study. Actual results may vary.





Cordless Power Flosser 3000

For those who want the very best for their mouth

Enjoy an effortless, thorough clean between your teeth and improve gum health. Our Quad Stream technology helps you cover more area with less effort, while Pulse Wave technology gently guides you from tooth to tooth




250ml Reservoir

The 250ml reservoir holds enough water for 60 seconds of cleaning. To refill either twist and remove, or use the side door to one handed refilling


Covers more area, effortlessly

X-shaped water streams cover more surface area for a faster, deeper clean versus manual floss



F3 Quad Stream nozzle

For maximum coverage and an extra thorough clean, use the Quad Stream nozzle to concentrate water into four wide streams between teeth and along the gum line. The softer rubber guide is gentle on gums and every nozzle simply clicks on and off



F1 Standard nozzle

Use the standard nozzle for a regular clean. It uses a single stream of water to remove debris and clean thoroughly between your teeth

Choose the perfect brush head for you



Each brush head fits every Philips Sonicare toothbrush



Sonicare For Kids
Mini



Sonicare For Kids
Tall



C1
ProResults



W2
Optimal White
(Black and White)



S2
Sensitive



G3
Premium Gum Care
(Black and White)



C3
Premium Plaque Defence
(Black and White)









A3
Premium All-in-One White
(Black and White)



All our brushes have in-built brushsync, when paired with a head showing this logo they will remind you when it is time to change your brush head

Which Sonicare is right for you



| | Sonicare 7100 series | Sonicare 6100 series | Sonicare 5300 series | Sonicare 3100 series | Sonicare For Kids | Cordless Power Flosser |
|---------------------------------|---|---|---|---|---|--|
| |  |  |  |  |  |  |
| Plaque removal | Up to 10x more vs. a manual toothbrush | Up to 10x more vs. a manual toothbrush | 7x more vs. a manual toothbrush | Gentle Clean, even for sensitive teeth and gums | Up to 75% more vs. a children's manual toothbrush ¹ | Up to 2x as effective as interdental brushes ² |
| Cleaning modes | 4 | 2 | 1 | 1 | 2 | 2 |
| Intensity settings | 3 | 3 | 2 | 1 | — | 3 |
| Sonicare App | ✓ | — | — | — | ✓ | — |
| 2-minute timer | ✓ | ✓ | ✓ | ✓ | ✓ | — |
| Pressure sensor | ✓ | ✓ | ✓ | ✓ | — | — |
| Accessories and Charging | Charging case, USB charger | Travel case, USB charger | USB charger | USB charger | — | Travel bag, USB-C charger |

Did you know



The best way to brush
is at a 45° angle
along the gum line



Slide and glide
with minimal
pressure





Scan me with your smart phone camera

www.philips.co.uk/sonicare