

PHILIPS



Oral Irrigator

Cordless Power Flosser 3000

Quad Stream technology

Pulse waves guide as you go



HX3828/22

Flossing reinvented with QuadStream technology

Removes upto 99% of plaque from below the gumline*

Effortlessly thorough cleaning between teeth to improve your gum health. Unique Quad Stream tip cleans more area with less effort for more effective flossing every time, so it's easy to get flossing right.

Easy and convenient

- Easy-to-fill 250 ml reservoir
- Effortless technique
- Standard tip for spot cleaning
- Comfort tip for sensitive gums

Faster, more effective flossing

- 2 flossing modes, 3 intensities
- Complete clean in 60 seconds
- Up to 150% more effective than floss in improving gum health
- Quad Stream technology for faster, more effective flossing

An effortless clean with advanced technology

- Quadstream technology: removes plaque from below the gumline

Cordless Power Flosser 3000

Oral Irrigator

HX3828/22

Highlights

Quad Stream technology




Unique X-shaped Quad Stream tip separates flow into 4 water jets, which cover more area between teeth and along the gum line. Get a faster, deeper clean without the hassle of string floss.

Cleans from below the gumline



The Quadstream technology removes up to 99% of plaque from 6mm deep periodontal pockets.

Up to 150% more effective than floss




Improves gum health in just 2 weeks****. and is upto 150% more effective than traditional floss***

2 modes, 3 intensities



Choose the clean that suits your smile. Clean mode delivers continuous flow for effective everyday cleaning. Deep Clean pulses for more thorough cleaning. Adjustable intensity level adapts for comfort.

Complete clean in 60 seconds




Access those hard-to-reach areas with a nozzle that rotates 360 degrees, so that you can use the device in any orientation. Get a thorough 360° clean in just 60 seconds from start to finish.

Comfort tip for sensitive gums




Want to take care of sensitive gums? Use the Comfort tip and gently guide the soft rubber tip along your gumline.

250ml reservoir




Reservoir capacity of 250 ml holds enough water for a recommended 60-second clean. Simply twist and remove or use the side-fill door for easy refilling.

Effortless technique



Just point, press and clean for more effective flossing. Aim the tip between your teeth and along the gumline — and let it do all the work for you!

Standard tip for spot cleaning



A standard single-stream tip boosts pressure for spot cleaning and removing food debris when you need it.

Cordless Power Flosser 3000

Oral Irrigator

HX3828/22

Specifications

Power

Voltage: Multi-voltage charger

Technical specifications

Energy Consumption: Off mode <0.3 W, Standby without display <0.5 W

Design and finishing

Color: White

Service

Warranty: 2-year limited warranty

Ease of use

Handle: Slim ergonomic design

Nozzle attachment: Easily clicks on and off

Performance

Cleaning: Clean entire mouth in 60-90sec

Items included

F1 Standard nozzle: 2

F3 Quad Stream nozzle: 2

Travel bag: 2

USB charging cable: 2

F2 Comfort nozzle: 2

USB wall adaptor: 1

Sonicare Cordless Power Flosser: 2

Nozzle case: 2

Modes

Clean: For a thorough everyday clean

Deep Clean: For a deep interdental clean

Intensities: 3 levels

Health benefits

Sub-gingival: Removes up to 99% of plaque from below the gumline**



* In an in-vitro study. Actual results may vary.

* *In a lab study (actual results may vary) with Quad Stream nozzle (in clean mode, setting 3), in the deepest 2mm of a 6mm periodontal pocket model.

* **between teeth, after two weeks, in patients with moderate to severe gingivitis, when used with Quad Stream nozzle on high setting with a manual toothbrush vs. using manual toothbrush with string floss

* ***when used with Quad Stream nozzle on high setting with a manual toothbrush in patients with moderate to severe gingivitis