PHILIPS
sonicare

1100

300015933351 (2025-08-26)

English

We hereby inform you that the user guides for Philips Sonicare electric toothbrush kits 1100 Series contain instructions for use of the models (types) of the electric toothbrushes listed in the table, which also includes their accessories, basic parameters, characteristics and the countries of production.

Model of Philips Sonicare kit consisting of electric toothbrush or toothbrushes with accessories	HX3901/01	HX3901/03	HX3902/01
Series	1100	1100	1100
Type of Philips Sonicare electric toothbrush (battery type, country of production, safety class)	HX390W1 (Li-Ion, made in China, Class III Constructions)	HX390MB (Li-Ion, made in China, Class III Constructions)	HX390W1 (Li-Ion, made in China, Class III Constructions)
Brushing modes	Clean	Clean	Clean
Brushing intensity levels	1	1	1
Series of Philips Sonicare toothbrush head (bristle stiffness, pc.)	i, InterCare (medium soft, 1 pc.)	i, InterCare (medium soft, 1 pc.)	i, InterCare (medium soft, 1 pc.)
USB-A cable	1 pc.	1 pc.	1 pc.
Toothbrush head reminder mode	No	No	No

Storage and Usage Conditions

Transport and storage temperature: -20 °C to +60 °C (-4 °F to 140 °F)

Transport and storage relative humidity: 10% - 90%

Temperature of usage: 0 °C - 40 °C (32 °F to 104 °F)

Relative humidity of usage: 20% - 90%

The manufacturer: Philips Consumer Lifestyle B. V.
Tussendiepen 4, 9206 AD, Drachten, The Netherlands
For household needs

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