

PHILIPS

sonicare



Rechargeable
toothbrush

Series 6100

10 x more plaque removal*

Gentle on sensitive teeth and
gums

Next-Generation Sonicare
Technology

Visual Pressure Sensor



HX7401/11

Advanced care for sensitive teeth and gums

Gently removes 10 x more plaque*

Thorough cleaning can still be gentle on sensitive teeth and gums with 10 x more plaque removal*. Featuring our Next-Generation Sonicare Technology for a powerful yet gentle clean, even in those hard-to-reach areas.

Brushing made effortless

- Choose your ideal cleaning experience
- Guided brushing sessions
- Keep brushing your best with the brush head replacement reminder

Designed for your lifestyle

- 21 days of regular brushing

Feel the care of Sonicare

- Extra gentle on sensitive teeth
- Next-Generation Sonicare Technology
- Help protect your gums with the pressure sensor

Series 6100

Rechargeable toothbrush

HX7401/11


Highlights

Extra care for sensitive teeth




Thorough cleaning can still be gentle on sensitive teeth and gums. This S2 Sensitive brush head features a unique design with long, thin, ultra-soft bristles that keep every brush stroke kind. With more than 3000 densely packed bristles, it helps you remove 10 x more plaque than a manual toothbrush*.

Next-Gen Sonicare Technology



This electric toothbrush uses our Next-Generation Sonicare Technology to give you a truly even clean throughout your whole mouth. The motor adjusts power automatically so there is a zero drop in performance when you reach those harder-to-clean areas. Enjoy a gentle yet effective clean for your teeth and gums with up to 62,000 bristle movements. Sonicare Fluid Action supports the bristles to clean by driving fluid deep between teeth and along the gumline.

Visual pressure sensor




This Sonicare electric toothbrush has a light ring at its base that lets you know if you're overdoing the pressure. Simply ease off when it lights up and your gums will stay protected.

Modes and Intensities




Give your clean a boost with 6 brush settings. Whether you want a little more vigour or a specific focus for your clean, you can have it. Choose between Clean and Sensitive mode. Pick one of three intensity settings.

SmartTimer and BrushPacer



Ace your session time with Sonicare's cleaning timers. Every 20 seconds, the BrushPacer will prompt you to clean a new area. After 2 minutes, the SmartTimer will indicate your session is complete.

Replacement reminder



Did you know brush heads become less effective after three months? Dental professionals recommend changing your brush head regularly. That's why Philips Sonicare toothbrushes have brush head replacement reminder. It tracks how often and how hard you brush, then reminds you to replace your brush head once it's time for a new one.

21-day battery life



Experience up to 21 days of regular brushing after a single charge, bringing a new level of convenience to your brushing routine.

Series 6100

Rechargeable toothbrush

HX7401/11

Specifications

Intensities

High: To boost your clean
Medium: For an everyday clean
Low: For sensitive teeth and gums

Power

Voltage: 100-240 V

Technical specifications

Battery: Rechargeable
Operating time (full to empty): 21 days
Battery type: Lithium ION
Energy consumption: Standby without display
0.08 W

Design and finishing

Colour: Black

Service

Warranty: 2-year limited warranty

Ease of use

Handle: Slim ergonomic design
Battery indicator: Illuminated icon indicates battery life
Handle compatibility: Easy click-on brush heads
Timer: BrushPacer and Smartimer

Items included

Handle: 1 6100 rechargeable toothbrush
Brush head: 1 S2 Sensitive
Charger: 1 USB-A charger

Cleaning performance

Plaque removal: Up to 10 x more effective*
Speed: 62,000 brush movements/min

Modes

Clean: For exceptional everyday cleaning
Sensitive: For sensitive teeth and gums

Smart sensor technology

Pressure feedback: Light ring illuminates purple, Vibration and pulsing sound
Replacement Reminder: BrushSync lets you know when to replace your brush head

