## EasyStart

This Philips Sonicare model comes with the EasyStart feature activated. The EasyStart feature gently increases the power over the first 14 brushings to help you get used to brushing with the Philips Sonicare. The EasyStart feature works with Clean mode only. To deactivate EasyStart see "Activating or deactivating features".

## Smartimer

The Smartimer indicates that your brushing cycle is complete by automatically turning off the toothbrush at the end of the brushing cycle. Dental professionals recommend brushing no less than 2 minutes twice a day.

## QuadPacer

The Quadpacer is an interval timer that beeps every 30 seconds to remind you to brush all teeth evenly by dividing your mouth into 4 sections (Fig. 7).

## Pressure sensor

Your Sonicare toothbrush measures the pressure you apply while brushing to protect your gums and teeth from damage.

If you apply excess pressure, the handle changes its vibration until you reduce the pressure.

## Activating or deactivating features

You can activate or deactivate the following features of your toothbrush:

- EasyStart
- Pressure Sensor Feedback

To activate or deactivate these features, follow the instructions below:

## EasyStart

- Put the handle on the plugged-in charger.
- Keep the power button pressed until you hear a single short beep (after 3 seconds).
- Release the power button.
  - Triple tone of low-medium-high means the EasyStart feature has been activated. The battery indicator blinking green twice shows the EasyStart is activated.
  - Triple tone of high-medium-low means the EasyStart feature has been deactivated. The battery indicator blinking amber once shows the EasyStart is deactivated.

**Note:** To achieve optimal efficacy, EasyStart needs to be deactivated.

## Pressure Sensor Feedback

- Put the handle on the plugged-in charger.
- Keep the power button pressed until you hear a series of three short beeps and see the Pressure Sensor icon light up (after 7 seconds) .
- Release the power button.
  - Triple tone of low-med-high means the Pressure Sensor Feedback feature has been activated. The battery indicator light will blink green twice to confirm activation.
  - Triple tone of high-med-low means the Pressure Sensor Feedback feature has been deactivated. The battery indicator light will blink amber once to confirm deactivation.

## Cleaning

You should clean your Philips Sonicare toothbrush regularly to remove toothpaste and other residue. Failure to clean your Philips Sonicare toothbrush may result in an unhygienic product and damage to it may occur.

**Warning: Do not clean brush heads, the handle or the charger in the dishwasher.**

**Warning: Do not clean with boiling hot water.**

## Toothbrush handle

- Remove the brush head and rinse the metal shaft with warm water. Make sure you remove any residual toothpaste on the handle and around the buttons at least once a week, but preferably after each use (Fig. 9).

**Caution: Do not push on the rubber seal on the metal shaft or around the buttons with sharp objects, as this may cause damage.**

- Wipe the entire surface of the handle with a damp cloth.

**Caution: Do not tap the handle on the sink or counter to remove excess water.**

- Dry the handle with a soft cloth or towel.

## Brush head

- Rinse the brush head and bristles after each use and allow to air dry (Fig. 10).
- Remove the brush head from the handle and rinse the brush head connection with warm water at least once a week.

## Charger

- Unplug the charger before you clean it.
- Wipe the surface of the charger with a damp cloth.
- Dry the charger with a soft cloth or towel.

## Storage

If you are not going to use the product for a long time, unplug it, clean it (see chapter "Cleaning") and store it in a cool and dry place away from direct sunlight.

Charge the toothbrush regularly to prevent the battery from discharging completely. A completely discharged battery may fail to charge again.

## Brush head replacement

Replace Philips Sonicare brush heads every 3 months to achieve optimal results. Use only Philips Sonicare replacement brush heads.

- Replace the brush head if signs of wear appear. Stop using a brush head with crushed or bent bristles.
- Smart fading bristles provide a visible replacement reminder: when the blue bristles fade to white at one-third of their length, replace the brush head.

## Locating the model number

Look on the bottom of the Philips Sonicare toothbrush handle for the model number HX507x.

## Troubleshooting

This chapter summarizes the most common problems you could encounter with the appliance. If you are unable to solve the problem with the information below, visit **www.philips.com/support** for a list of frequently asked questions or contact the Consumer Care Center in your country.

Problem	Possible Cause	Solution
I cannot fully attach the brush head. There is a gap between the brush head and handle.		This gap is normal and is needed for the brush head to vibrate appropriately. The brush head must be able to move to create the right amount of vibrations.
My Philips Sonicare toothbrush vibration is less powerful than before.	The brush head is too close to the body of the toothbrush or the toothbrush needs to be charged.	Remove the brush head from the handle and re-set it on the metal shaft leaving a small gap. (1-2mm). <p>When charging, make sure the charger is plugged into a power supply and placed on a flat surface before placing the toothbrush on the charger. Allow the toothbrush to charge for 10 hours.</p>
I charged my toothbrush, but it only works for a short time.	The base of the toothbrush is not placed flat on the charger when attempting to charge.	Make sure the charger is plugged into a power supply and placed on a flat surface before placing the toothbrush on the charger. Allow the toothbrush to charge for 10 hours.

## Warranty and support

Philips offers a 2-year warranty. Warranty exceptions apply. This warranty does not affect your rights as a consumer under law. For more information, visit **www.philips.com/support** or call us.

## Warranty restrictions

The terms of the international warranty do not cover the following:

- Brush heads.
- Damage caused by use of unauthorized replacement parts.
- Damage caused by misuse, abuse, neglect, alterations or unauthorized repair.
- Normal wear and tear, including chips, scratches, abrasions, discoloration or fading.

## Recycling

- This symbol means that electrical products and batteries shall not be disposed of with normal household waste (Fig. 11).
- Follow your country's rules for the separate collection of electrical products and batteries.

## Removal of built-in rechargeable battery

CONTAINS LITHIUM ION RECHARGEABLE CELL. MUST BE RECYCLED OR DISPOSED OF PROPERLY.

The built-in rechargeable battery must only be removed by a qualified professional when the appliance is discarded. Before removing the battery, make sure that the appliance is disconnected from the wall socket and that the battery is completely empty.

**Take any necessary safety precautions when you handle tools to open the appliance and when you dispose of the rechargeable battery. Be sure to protect your eyes, hands, fingers, and the surface on which you work.**

**If batteries are damaged or leaking, avoid contact with the skin or eyes. If this occurs, immediately rinse well with water and seek medical care.**

**When you handle batteries, make sure that your hands, the product and the batteries are dry.**

**To avoid accidental short-circuiting of batteries after removal, do not let battery terminals come into contact with metal objects (e.g. coins, hairpins, rings). Do not wrap batteries in aluminum foil. Tape battery terminals or put batteries in a plastic bag before you discard them.**

Only remove the rechargeable battery when you discard the product. Before you remove the battery, make sure that the product is disconnected from the wall socket and that the battery is completely empty.

- Remove the inner body from the handle.
- Remove the rechargeable battery from the inner body.

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