Neonatal and Infant Single-Patient NBP Cuffs

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

EN - Neonatal and Infant Single-Patient NBP Cuffs Instructions for Use

Associated Documents - For more information, refer to the *Air Hoses for Noninvasive Blood Pressure Cuffs Instructions for Use*. **Intended Use** - The Philips Neonatal and Infant Single-Patient NBP Cuffs are to be used with compatible neonatal monitors (see *Compatibility*) intended for use by, or under the supervision of, a licensed physician or other healthcare provider for the noninvasive measurement of neonatal and infant human blood pressure. **Intended Purpose** - In accordance with European Medical Device Regulation 2017/745, this accessory device enables a vital signs patient monitoring instrument to record and generate alarms for multiple physiological parameters when monitoring neonates. **Indications** - Philips Neonatal and Infant Single-Patient NBP Cuffs are indicated for the measurement of NBP and pulse rate for monitoring purposes. These devices are limited by the indications for use of the connected monitor in healthcare facilities. These devices are intended to interact with patient intact skin only. **Contraindications** - There are no known contraindications. **Compatibility** - Neonatal NBP Cuffs can be used with any monitors that list the cuffs as accessories in that product's instructions for use. The Neonatal NBP Cuffs can be used with the following Neonatal Blood Pressure Air Hoses: 989803209741 or M1596C (both 1.5 m Hoses) and 989803209751 or M1597C (both 3.0 m Hoses).

Before use, read all warnings and cautions listed in this IFU. Also refer to the IFU for the monitor or defibrillator for additional warnings and cautions.

Warnings - Single-patient cuffs are not intended for use on more than one patient. • Use only monitor hoses identified for use with neonatal cuffs. Use of other hoses can prevent the monitor from obtaining readings. • Be sure to select the correct patient type setting on the monitor. Do not apply the higher adult inflation, overpressure limits, and measurement duration to neonatal patients. This may cause a patient injury or inaccurate measurements. • Be sure that the cuff is correctly matched to patient type and size (limb circumference ranges are printed on each cuff). Cuffs that are too small may produce false high blood pressure readings; cuffs that are too large may produce false low-pressure readings. • Kinking of the tube can result in inaccurate readings: avoid placement of the cuff and tubing which could lead to kinking. Kinking may cause patient injury. • To prevent damage to the cuff (for example, leaks), shut off auto function or remove the cuff from monitor when not on patient. • Use cuff only under direct supervision by a trained healthcare professional when attached to automated monitors without fail-to-read alarms. • Periodically inspect cuff, connectors, and hose for signs of wear or deterioration. Replace damaged or deteriorated components, as required. • Replace cuff if hook and loop fastener fails to hold during inflation. • Follow your institution's guidelines for long-term patient care. • Do not sterilize cuffs. Do not clean or disinfect. Do not use bleach. • Modification of this equipment is not allowed.

Cautions - Federal (US) law restricts this device to sale by or on the order of a medical practitioner. • Do not substitute any component of the NBP product since a substitution might result in a measurement error.

Applying the Single Patient NBP Cuff - Notes: • Whenever visual inspection reveals a cuff is no longer suitable for continued use, does not perform as intended, or severe soiling occurs, dispose of the cuff. Refer to *Disposal* section. • Refer to the figures on page i. 1-Select the proper cuff size for the patient's limb circumference. Refer to *Select Appropriate Cuff Size*. 2-Apply the cuff to the limb such that the side displaying the Δ is toward the patient. The arrow labeled "Arteria" must be over the limb artery. 3-Wrap the cuff snugly around the patient's limb. Ensure that the index edge falls within the range indicated by the range arrow. 4-To prevent leaking, secure the connector(s) on the cuff to the air hose connector. **Select Appropriate Cuff Size** - NBP Cuffs are available to suit a variety of patient applications and are color-coded for easy size selection:

	Neonatal # 1	Neonatal # 2	Neonatal # 3	Neonatal # 4	Neonatal # 5
Single-Patient Cuffs	989803167181, M1866B 989803167231, M1866S	989803167191, M1868B 989803167241, M1868S	989803167201, M1870B 989803167251, M1870S	989803167211, M1872B 989803167261, M1872S	989803167221, M1873B 989803167271, M1873S
Circumference (cm)	3.1–5.7	4.3–8.0	5.8–10.9	7.1–13.1	10.0–15.0
Cuff Text Color	Orange	Light Blue	Green	Navy	Burgundy

Reordering Information - The following products are CE marked and available in the European Union: 989803167181 (M1866B) — Neonatal Single-Patient Cuff Size #1; 989803167231 (M1866S) — Neonatal Soft Single-Patient Cuff Size #1; 989803167191 (M1868B) — Neonatal Single-Patient Cuff Size #2; 989803167241 (M1868S) — Neonatal Soft Single-Patient Cuff Size #2; 989803167201 (M1870B) — Neonatal Single-Patient Cuff Size #3; 989803167251 (M1870S) — Neonatal Soft Single-Patient Cuff Size #3; 989803167211 (M1872B) — Neonatal Single-Patient Cuff Size #4; 989803167261 (M1872S) — Neonatal Soft Single-Patient Cuff Size #4; 989803167221 (M1873B) — Infant Single-Patient Cuff Size #5; 989803167271 (M1873S) — Infant Soft Single-Patient Cuff Size #5. **Cleaning and Care** - Neonatal NBP Cuffs are for single-patient use and are not intended to be cleaned or disinfected. If the cuff is soiled, discard it. Refer to the *Disposal* section that follows. **Disposal** - When finished using this single use cuff, follow approved disposal methods for medical waste as specified by your facility (hospital, clinic) or local regulations. **EU Incident Reporting** - Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established. **Environmental Specifications - Storage**: Temperature: -20° to +60° C (-4° to 140° F); Humidity: 0% to 90% RH@ 60° C (140° F); **Operating**: Temperature: 0° to +50° C (32° to 122° F); Humidity: 0 to 95% RH.

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Copyright ©2022
Böblingen GmbH