

Philips Avent single electric breast pump

Philips Avent is dedicated to helping mothers breastfeed for longer, because it is proven that breast milk is vitally important for healthy infant development

Research based product development. We have collaborated with leading scientists and healthcare professionals for over 25 years and conducted research with more than 15,000 mothers worldwide.

What makes this breast pump unique? Recent studies have proven a clear correlation between comfort and milk flow during expressing. Relaxed mothers feel more comfortable and express significantly more milk more easily. Independant research also confirms that infant suckling is a synergistic combination of breast compression and intra-oral vacuum, originating from the peristaltic action of the infant's tongue. The Philips Avent Comfort breast pump is based on exactly these insights.





Unique, natural expressing position

The breast pump has a unique breast shield design. Mothers can sit in a comfortable position when expressing, with no need to lean forward yet milk still flows directly into the container. This improved position helps mothers feel more comfortable and, as a result, express more milk, naturally.

93% of mothers (strongly) agreed that being comfortable when expressing is important for good milk flow.*



Clinically proven to be effective

The breast pump mimics the infant's natural peristaltic suckling behaviour by combining breast compression with gentle suction. It features our clinically proven massage cushion with soft petals that gently compress the breast to stimulate milk flow. This innovative cushion has a distinctive, velvety texture that feels warm on the breast, further helping mothers to relax whilst expressing. The removable cushion also prevents milk pooling, aiding hygiene and eliminating the risk of lost milk.

92% of mothers were (very) satisfied with the effectiveness of the breast pump.*



Simple settings

When switched on, the pump automatically starts in gentle stimulation mode to help trigger let down. Then the mother can choose between 3 pumping settings to make milk flow most comfortable and effective for her. Easy to use on the go with batteries.

98% of mothers (strongly) agreed that the pump settings are easy to use.*

*Independent home placement test among 85 mothers in the UK, August 2012

Product components

Comfort breast pump	Specifications	Benefits
Pumping mode	Single electric	Effortless expressing at the touch of a button
Key features	Unique compact breast shield design with integrated funnel Clinically proven massage cushion with distinctive, velvety texture Simple settings, usable with batteries	Unique, natural expressing position for easier milk flow Gentle stimulation of natural let down and milk flow Stimulation mode plus 3 expression settings to choose from. Easy to use on the go
Vacuum levels	Stimulation mode: 128mmHg. Expression setting 1: 169mmHg Expression setting 2: 209mmHg. Expression setting 3: 250mmHg	Mothers can choose their most comfortable expressing setting
What's in the pack	1x breast pump expression unit (including 1x Philips Avent Natural bottle) 1x motor unit including tubing & connecting cap 1x newborn flow teat 1x hygiene cover for cushion 1x sealing disc for breast milk storage 1x power adapter 1x spare diaphragm 1x instruction manual for Comfort breast pump 1x instruction manual for Natural bottle 1x day breastpad sample pack 1x night breastpad sample pack	All that mothers need to start expressing, feeding and storing breast milk
Materials	Polypropylene: pump body, hygiene cover for funnel, Natural bottle. Silicone: cushion, diaphragm, valve, newborn flow teat TPE: sealing disc Entirely BPA free	Safe and hygienic to use. Milk can never get into the tubing or motor unit. All parts can be sterilised except for the electrical parts.
Dimensions	Expression unit: Total height: 17.1cm (when fully assembled) Expression unit: Total width: 13.5cm (when fully assembled) Motor unit: Total height: 6.8cm, Total width: 14.2cm, Total depth: 10.4cm	Compact size for easy expressing. Reduced storage space.
Weight	Expression unit: 125g (when fully assembled) Motor unit including tubing: 307g	Lightweight for easy transport
Cleaning	All parts can be immersed in water and sterilised except for the electrical parts	Easy to keep the pump hygienic. Breast milk can never get into the tubing or motor unit.
Assembly	Easy visual matching of parts. Few small pieces	Intuitive assembly
Instructions of use	Printed instruction manual with easy-reference troubleshooting guide	All a mother needs to know about starting to express in one leaflet
Warranty	2 years	Philips Avent consumer care service available for any support
Accessories	Large soft massage cushion available for purchase separately. Funnel size: Ø2.5cm (Ø1.95cm funnel size on regular cushion that comes standard with the Philips Avent Comfort breast pump) Avent Comfort breast pump)	Comfortable expressing for larger than average nipples
Compatibility	Compatibility across the Philips Avent range (including e.g. Classic bottles and milk storage containers)e.g. Classic bottles and milk storage containers)	Flexibility in combining different feeding and storage options
Available spare parts	Pump body, massage cushion, massage cushion large, diaphragm and valve, tubing, power adapter, motor unit, hygiene cover for cushion, battery lid on motor unit	Can be ordered online or via the Philips Avent consumer care service

Other products to help mothers breastfeed for longer



2 in 1 thermopads: can be used warm, to stimulate milk flow, or cold to soothe sore and engorged breasts.



Nipple protectors: ultra-fine, tasteless and odourless. Allows mothers to keep breastfeeding in case of sore or cracked nipples. Available in small and standard size.



Day and night breast pads: ultra dry, with silky soft feel. Designed to be extra discreet keeping mothers comfortable between feeds. Available for day and night.



Milk storage containers: for easy and safe storage of breast milk in the fridge or freezer. Compatible with all Philips Avent breast pumps.

Moistu sensiti delicat

Moisturising nipple cream: softens dry or sensitive nipples and helps prepare the delicate breast area for feeding.



