

Patient interface selection



Ventilator compatibility

Single limb circuit	EE elbows (leak 1, 2)	•	EE elbow	EE elbow	EE elbow	EE elbow	•	•	•	•
Dual limb circuit	SE elbow		SE elbow	SE elbow	SE elbow	SE elbow	•	•		
Patient treatments										
Bronchoscopy procedure	Bronchoscopy elbow		Bronchoscopy elbow							
Aerosol medication delivery			NIVO elbows							
Patient characteristics										
Claustrophobia	•	•								•
Facial irregularities	•	•								•
With no teeth	•	•	•	•	•	•	•	•	•	•
Easier mouth access			with CapStrap	•	with CapStrap		•	•	•	•
Mouth breathing	•	•	•	•	•	•				
Application indications										
Rapid application	•	•	with CapStrap	•	with CapStrap			•		
Frequent NIV patient	•		•	•	•	•	•	•		•
Post-surgical/recovery			•	•	•	•				
Headgear type										
Four-point	•	•	•	•	•	•	•			
Four-point with chin strap										•
CapStrap			•	•	•			•		
Three-point						•				
Bonnet									•	

Elbows

