




## Phoenix guidewire

### Technical specifications



Part number	Type	Diameter	Length	Tip	Style
PG14300LF	Silicone coated nitinol core	.014"	300 cm	Floppy	Light support
PG14300XF	Silicone coated stainless steel core	.014"	300 cm	Floppy	Extra support

### Indications

The Phoenix guidewire family is designed for percutaneous entry into the peripheral vascular system as an intravascular device. It gently navigates the vascular system while maintaining enough strength and rigidity to enable a catheter to be directed over the guidewire with ease.

### Contraindications

This device is not intended for use in the coronary or neurovascular arteries.

***The Phoenix guidewire is compatible with the Phoenix atherectomy system. This product is currently only available in the United States of America.***