

PHILIPS
RESPIRONICS

Sleep and Respiratory Care

PAP battery kit with
integrated uninterruptible
power supply (UPS)

Respiratory battery kit



Give your patients the **confidence and convenience** of powering their Philips Respironics PAP or respiratory therapy device with our lithium ion battery kits.

All DreamStation and System One PAP devices, as well as Trilogy ventilators, can be powered with the Philips Respironics lithium ion battery kits.

Peace of mind

Your patients will know that when the power goes out, the battery pack's integrated uninterruptible power supply (UPS) will help reduce therapy disruptions.

Traveling convenience

It's our lightest and smallest battery pack ever, offering easy portability so patients can continue therapy away from home. It's great for camping, hotel use, or active duty military. The PAP battery kit uses the power supply from the patient's Philips Respironics sleep therapy device to reduce the number of cords needed for travel.

Confidence for extended trips

It meets FAA requirements for airplane travel (90W lithium ion battery pack), and it provides more than 14 hours of battery life.*

Capabilities for PAP patients

Continuously connecting the PAP battery kit to the patient's Philips Respironics sleep therapy device's power supply keeps it charged and ready for travel or UPS protection.

Capabilities for respiratory patients

- The battery pack provides more than 8 hours of battery life.**
- It easily attaches to a wheelchair via the carrying case and its attachment strap.

*PAP run-time test conditions: CPAP mode, pressure 10 cm H₂O, no heated humidification or heated tubing.

**Respiratory run-time test conditions for Trilogy: S/T mode, IPAP 15 cm H₂O, EPAP 5 cm H₂O, BPM 12, inspiratory time 0.8 s, passive, compliance 20, resistance 20.

Benefits for patients

- The battery pack provides more than 14 hours of battery life for PAP devices*.
- It's compact and lightweight at just 2.25 lbs.
- It meets FAA requirements for airplane travel.
- It comes with a six-month warranty.
- Its lithium ion battery is environmentally conscious.



No prescription is required.



Ordering information for PAP

Description	Part number
PAP Lithium Ion Battery Kit*	1124009
Includes a rechargeable lithium ion battery pack, a PAP device cable, 2 adapters (for use with System One 50 series only), and a battery case.	
Supports DreamStation, System One 60 series, System One 50 series CPAPs and BiPAPs. Utilizes DreamStation or System One power supply. Additional power supply not included.	
Replacement parts	Part number
Li Ion battery, PAP device cable	1121736
Li Ion Battery, Cable Adapter (50 series adapter, qty 2)	1127799
Li Ion Battery Case	1126947
DreamStation 65W power supply	1120136
DreamStation 80W power supply	1118499
DreamStation/System One power cord	1005894
System One 60W power supply (for H and HS models)	1091398
System One 80W power supply (for TS model only)	1091399



Ordering information for Respiratory

Description	Part number
Respiratory Lithium Ion Battery Kit**	1124010
Includes a rechargeable lithium ion battery pack, a Trilogy cable, a 65W power supply, a battery case, and a battery case attachment strap.	
Supports Trilogy as a battery backup. UPS feature not available for use in this kit.	
Replacement parts	Part number
Li Ion battery, Trilogy cable	1121737
Li Ion battery, BiPAP (A30/A40) device cable (Canada only)	1121879
Li Ion Battery Case	1126947
Li Ion Battery Case Attachment Strap	1126941



Lithium ion battery pack specifications

Warranty	Six-month
Dimensions	1.25" (3.175 cm) x 5.25" (13.335 cm) x 9.0 (22.86 cm)
Weight	2.25 lbs. / 1.02 kg
Prescription	Not required
Capacity	90W
Charge/discharge indication	Push-button illumination
Uninterruptible power supply (UPS)	Yes — for use with PAP lithium ion battery kit only

Battery technology	Lithium ion
Input voltage range	11.0-14.4 VDC
Output voltage range	PAP use: 11.0-15.0 VDC Respiratory use: 10.5-12.2 VDC
Run time	PAP use: 14 hours*** Respiratory use: 8 hours****
Minimum life cycle	≥70% of rated capacity after 1000 cycles
Recharge time	3-4 hours

For more information, call 1-800-345-6443 or visit www.philips.com/sleeptherapy.

*PAP Lithium Ion Battery Kit for use with the Philips Respironics DreamStation and System One CPAPs, BiPAPs, autoSV, BiPAP AVAPS, and BiPAP S/T sleep therapy devices.

**Respiratory Lithium Ion Battery Kit for use with Trilogy.

***PAP run-time test conditions: CPAP mode, pressure 10 cm H₂O, no heated humidification or heated tubing.

****Respiratory run-time test conditions for Trilogy: S/T mode, IPAP 15 cm H₂O, EPAP 5 cm H₂O, BPM 12, inspiratory time 0.8 s, passive, compliance 20, resistance 20.

© 2016 Koninklijke Philips N.V. All rights reserved.
Specifications are subject to change without notice.

www.philips.com/respironics



edoc 7/25/16 MCI 4107651
1010 Murry Ridge Lane, Murrysville, PA 15668
800-345-6443 • 724-387-4000