



The first choice for comfort and fit

Respironics PerforMax total face mask for multi-patient use

The Respironics PerforMax supports comfort, may help skin break-down, and eliminates nasal bridge challenges by sealing around the perimeter of the face, where patients have less pressure sensitivity and smoother facial contours. The PerforMax also simplifies mask selection and speeds setup by fitting a wide range of patients with only three sizes.

Key advantages

- Pediatric to adult extra large size fits most patients
- Perimeter facial contact supports comfort and simplifies fitting
- The clear, unobstructed view of the PerforMax may help to reduce patient claustrophobia

PHILIPS
sense and simplicity

Respironics PerforMax

The Respironics PerforMax for multi-patient use has additional advantages to support patient comfort.

- Interchangeable elbows accommodate quick transitions between single- and dual-limb ventilator circuits.
- Snap clips require minimal effort to attach and remove headgear.

Description	Quantity	Part number
PerforMax EE/SE, small*, multi-use	1 mask/4 headgear	1052553
PerforMax EE/SE, large, multi-use	1 mask/4 headgear	1052554
PerforMax EE/SE, extra large, multi-use	1 mask/4 headgear	1080420
PerforMax EE, small*, multi-use	1 mask/4 headgear	1052555
PerforMax EE, large, multi-use	1 mask/4 headgear	1052556
PerforMax EE, extra large, multi-use	1 mask/4 headgear	1080421
PerforMax SE, small*, multi-use	1 mask/4 headgear	1052568
PerforMax SE, large, multi-use	1 mask/4 headgear	1052569
PerforMax SE, extra large, multi-use	1 mask/4 headgear	1080422
PerforMax EE/SE, small*, multi-use	1 mask/1 headgear	1044385
PerforMax EE/SE, large, multi-use	1 mask/1 headgear	1044386
PerforMax EE/SE, extra large, multi-use	1 mask/1 headgear	1080423
PerforMax entrainment elbow (EE)**	1	1045313
PerforMax standard elbow (SE)**	1	1043210
PerforMax bronchoscopy elbow**	5	1072652
PerforMax bronchoscopy elbow**	1	1072651
PerforMax small headgear	5	1051857
PerforMax large headgear	5	1051816
PerforMax extra large headgear	5	1080518
PerforMax small headgear	9	1045325
PerforMax large headgear	9	1045314
PerforMax extra large headgear	9	1080519
Accessory swivel replacement kit	1	7041
PerforMax standalone sizing gauge (small, large)	1	1052529
Nasogastric tube sealing pad	10	452040
PerforMax cleaning and disinfection guide	1	1066982

* Intended for use by pediatric patients 7 years or older and > 20 kg

** Mask not included.



Leak 2 entrainment elbow (EE)



Standard elbow (SE)



Compatible with bronchoscopy elbow

Please visit www.philips.com/performax



© 2011 Koninklijke Philips Electronics N.V.
All rights are reserved. PerforMax and Respironics are trademarks of Respironics, Inc., dba Philips Respironics, and its affiliates.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Philips Healthcare is part of Royal Philips Electronics

www.philips.com/healthcare
healthcare@philips.com

Printed in The Netherlands
4522 962 76921 *AUG 2011