



Clear, reliable and fade-resistant



Philips medical recording paper

Protect your equipment investment, your records and your trace recording quality with Philips medical recording paper.

Designed and validated specifically for use with Philips patient monitors, defibrillators, and cardiographs, Philips medical recording paper is manufactured to exacting standards. The result is clear trace recordings that have exceptional archivability, without the premature fading you might experience with other papers.

In fact, recordings on our chemical/thermal paper for our fetal monitors, bedside monitors and defibrillators can last as long as five years, and Philips provides a 30 year warranty on our PTP brand anti-fade paper when used with specific models of Philips cardiographs.

In addition to providing excellent trace quality and fade-resistance, the unique chemical makeup of Philips medical

recording paper avoids build-up on print heads, which decreases the threat of damage, helps to prevent paper jams that can disrupt your workflow, and increases the time between print head cleaning. Philips paper also has special markings that make it easy to line up the paper.

We're so confident in the excellence of our paper that some Philips monitor, defibrillator, and cardiograph warranties are voided if print head damage is caused by use of a competitive paper. To receive the quality, reliability and performance you expect from Philips, use Philips medical recording paper.

PHILIPS
sense and simplicity

Fetal monitoring

Prod #	Description	Used with	Type	Dimensions	Packaging
92700484	Chemical/thermal paper - No header - Green grid	8030A 8031A 8040A 8041A 8232A/B	Z-fold	151 x 100mm Total length = 15m	1 pack = 150 shts 1 case = 40 packs 6,000 sheets total
92700485	PPR-ECG chart paper - No header - Orange grid, 1 min. (3cm highlights)	8030A 8031A 8040A 8041A 8232A/B	Z-fold	151 x 100mm Total length = 15m	1 pack = 150 shts 1 case = 40 packs 6,000 sheets total
92700630050	Chemical/thermal paper - No header - 2 channels - Green grid, KPA scale	8030A 8031A 8040A 8041A	Z-fold	151 x 100mm Total length = 15m	1 pack = 150 shts 1 case = 50 packs 7,500 sheets total
M1910A	Chemical/thermal paper - No header - Green grid, KPA scale	M1350A/B/C M1351A M1353A M2702A M2703A M2704A M2705A M2925A	Z-fold	150 x 100mm Total length = 15m	1 pack = 150 shts 1 case = 40 packs 6,000 sheets total
M1191A	Chemical/thermal paper - No header - Green grid, KPA scale	M1350A/B/C M1351A M1353A M2702A M2703A M2704A M2705A M2925A	Z-fold	150 x 100mm Total length = 15m	1 pack = 150 shts 1 case = 40 packs 6,000 sheets total
M1913A	Chemical/thermal paper - No header - Green grid	M1350A/B/C M1351A M1353A M2702A M2703A M2704A M2705A M2925A	Z-fold	150 x 100mm Total length = 15m	1 pack = 150 shts 1 case = 40 packs 6,000 sheets total
M1913J	Chemical/thermal paper - No header - Green grid for Japan	M1116A/B M3176B	Z-fold	150 x 100 mm Total length = 15m	1 pack = 150 shts 1 case = 40 packs 6,000 sheets total
M2206A	PTP-brand anti-fade, thermal paper for PageWriter XL - 2 7/8 in. header - Red grid	M1700A M1701A M1702A	Z-fold	8.5 x 11 in. Total length = 183 ft.	1 pack = 200 shts 1 case = 8 packs

Patient monitoring and recording modules

Prod #	Description	Used with	Type	Dimensions	Packaging
40469A	Chemical/thermal array paper - No header - Multi-channel - Unprinted	M1117A	Z-fold	120 x 150mm Total length = 45m	1 pack = 300 shts 1 case = 20 packs 6,000 sheets total
40477A	Chemical/thermal array paper - No header - Unprinted	M1116A/B M3176B	Roll	50mm Total length = 30m	1 box = 10 rolls
40477B	Chemical/thermal array paper - No header	M1116A/B M3176B	Roll	50mm Total length = 30m	1 box = 10 rolls 1 case = 8 boxes 80 boxes total
M4816A	Chemical/thermal array paper - No header	M3925A M3176C A1/C1/A3/C3	Roll	50mm Total length = 30m	1 box = 10 rolls
M4817A	Chemical/thermal array paper - No header	M3925A M3176C A1/C1/A3/C3	Roll	50mm Total length = 30m	1 box = 10 rolls 1 case = 8 boxes 80 boxes total
PSE11268	Seiko recording paper - No header - 4 channel	M3160A	Roll	112mm Total length = 30m	1 box = 24 rolls
989803136891	Thermal array paper - No header	VS1, VS3, VM4, VM6, VM8	Roll	5.8cm Total length = 30.5m	1 box = 5 rolls
989803138131	Thermal array paper - No header	862146	Roll	6.35cm Total length = 6m	1 box = 20 rolls

Defibrillators

Prod #	Description	Used with	Type	Dimensions	Packaging
40453A	Chemical/thermal paper - 1-channel - Gray grid	43100A 43110A 43120A 78173A/W 78660A/B 78661A/B 78670/71A	Roll	50mm Total length = 30.5m	1 box = 10 rolls
40457C	Chemical/thermal array paper - 1-channel - 40mm gray grid	78560A Codemaster M1722A/B M1723A/B M2475B M3500B M3500B M3535A M3536A/J M4535A M5500B	Roll	50mm Total length = 30.5m	1 box = 10 rolls
40457D	Chemical/thermal array paper - 1-channel - 40mm gray grid	78560A Codemaster M1722A/B M1723A/B M2475B M3500B M3500B M3535A M3536A/J M4535A M5500B	Roll	50mm Total length = 30m	1 box = 10 rolls 1 case = 8 boxes 80 boxes total
989803138171	Chemical/thermal array paper - 1-channel - 65mm red grid	M3535A M3536A/J MRx OptB04	Roll	50mm Total length = 30.5m	1 box = 10 rolls
989803138181	Chemical/thermal array paper - 1-channel - 65mm red grid	M3535A M3536A/J MRx OptB04	Roll	50mm Total length = 30.5m	1 box = 10 rolls 1 case = 8 boxes 80 boxes total

Cardiography and stress systems

Prod #	Description	Used with	Type	Dimensions	Packaging
M2206A	PTP-brand anti-fade, thermal paper for PageWriter XL - 2 ⁷ / ₈ in. header - Red grid	M1700A, M1701A, M1702A	Z-fold	8.5 x 11 in. Total length = 183 ft.	1 pack = 200 shts 1 case = 8 packs 1,600 sheets total
M2207A	PTP-brand anti-fade, thermal paper for PageWriter XL - 6.7cm header - Red grid	M1700A, M1701A, M1702A	Z-fold	216 x 297mm Total length = 59.4m	1 pack = 150 shts 1 case = 40 packs 6,000 sheets total
M2485A	PTP-brand anti-fade, thermal paper for PageWriter XL - 2 ⁷ / ₈ in. header - Red grid		Z-fold	8.5 x 11 in. Total length = 92 ft.	1 pack = 100 shts 1 case = 16 packs 1,600 sheets total
M2486A	PTP-brand anti-fade, thermal paper for PageWriter XL - 6.7cm header - Red grid		Z-fold	216 x 297mm Total length = 29.7m	1 pack = 100 shts 1 case = 16 packs 1,600 sheets total
M1707A	Thermal paper for PageWriter XL - 2 ⁷ / ₈ in. header - Red grid	M1700A M1701A M1702A	Z-fold	8.5 x 11 in. Total length = 183 ft.	1 pack = 200 shts 1 case = 8 packs 1,600 sheets total
M1708A	Thermal paper for PageWriter XL - No header - Red grid	M1700A M1701A M1702A	Z-fold	8.5 x 11 in. Total length = 183 ft.	1 pack = 200 shts 1 case = 8 packs 1,600 sheets total
M1709A	Thermal paper for PageWriter XL - 6.7cm header - Red grid	M1700A M1701A M1702A	Z-fold	216 x 297mm Total length = 59.4m	1 pack = 200 shts 1 case = 8 packs 1,600 sheets total
M1710A	Thermal paper for PageWriter XL - No header - Red grid	M1700A M1701A M1702A	Z-fold	216 x 297mm Total length = 59.4m	1 pack = 200 shts 1 case = 8 packs 1,600 sheets total
M1737A	ECG paper for laser printer - 2.4 in. header - Red grid	860281 M1730A 43610B	Sheets	8.5 x 11 in.	1 pack = 500 shts 1 case = 10 packs 5,000 sheets total
M1738A	ECG paper for laser printer - 60mm header - Red grid	860281 M1730A 43610B	Sheet	216 x 297mm Total length = 59.4m	1 pack = 500 shts 1 case = 10 packs 5,000 sheets total
M2481A	Thermal paper for PageWriter - 2 ⁷ / ₈ in. header - Red grid	M1700-02A, M1770-72A, 860284 Touch TC-series	Z-fold	8.5 x 11 in. Total length = 92 ft.	1 pack = 100 shts 1 case = 16 packs 1,600 sheets total
M2483A	Thermal paper for PageWriter - 6.7cm header - Red grid	M1700-02A, M1770-72A, 860284 Touch TC-series	Z-fold	216 x 297mm Total length = 59.4m	1 pack = 100 shts 1 case = 16 packs 1,600 sheets total

Cardiography and stress systems

Prod #	Description	Used with	Type	Dimensions	Packaging
M3707A	Thermal paper for PageWriter - 2 7/8 in. header - Red grid	860252	Z-fold	8.5 x 11 in. Total length = 92 ft.	1 pack = 100 shts 1 case = 16 packs 1,600 sheets total
		860286			
		860287			
		860288			
		860290			
M3708A	Thermal paper for PageWriter - 6.7cm header - Red grid	860286	Z-fold	216 x 297mm Total length = 29.7m	1 pack = 100 shts 1 case = 16 packs 1,600 sheets total
		860287			
		860288			
		860290			
		860297			
9898031 37011	Thermal paper for StressVue - No header - Red grid	860295	Z-fold	8.5 x 11 in. Total length = 229 ft.	1 pack = 250 shts 1 case = 12 packs 3,000 sheets total
9898031 43031	Thermal paper for StressVue - No header - Red grid	860295	Z-fold	216 x 297mm Total length = 74.25m	1 pack = 250 shts 1 case = 12 packs 3,000 sheets total

**Philips Healthcare is part of
Royal Philips Electronics**

How to reach us

www.philips.com/healthcare
healthcare@philips.com

Asia

+49 7031 463 2254

Europe, Middle East, Africa

+49 7031 463 2254

Latin America

+55 11 2125 0744

North America

+1 425 487 7000

800 285 5585 (toll free, US only)

Please visit www.philips.com/medicalsupplies or call your local representative



© 2012 Koninklijke Philips Electronics N.V.
All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Printed in The Netherlands.
4522 962 77731 * FEB 2012