



Enabling special head applications

Philips dStream T/R Head coil solution for Ingenia 1.5T / 3.0T

Transmit/receive coil, consisting of a base, sliding coil and head support, that provides excellent spectroscopy results due to its higher B_1 field. In addition, it allows imaging of patients with stereotactic frames. The open design reduces claustrophobia, while ensuring good homogeneity.

- High B_1 field for excellent spectroscopy
- Spacious design to accommodate stereotactic frames

PHILIPS
sense and simplicity

dS T/R Head coil solution

Coverage	30 cm
Max. number of channels	single transmit single receive
Main applications	Brain, head Spectroscopy Patients with stereotactic frame
Coil solution type	Transmit/Receive
Data sampling	DirectDigital
Coil connection	T/R interface



Please visit www.philips.com/Ingenia



© 2011 Koninklijke Philips Electronics N.V.
All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Philips Healthcare is part of Royal Philips Electronics

www.philips.com/healthcare
healthcare@philips.com

Printed in The Netherlands
4522 962 70021 * JAN 2011