



**PHILIPS**

FlexTrak Mammo

MR coils

## Designed to **improve your efficiency**

FlexTrak Mammo workflow solution for Ingenia 1.5T and 3.0T

### Key benefits

---

- Dedicated breast patient transport system
  - Available with 7ch and 16ch dS Breast coils
  - Workflow solution offering high quality imaging and patient comfort
- 

FlexTrak Mammo and the wide bore Ingenia MR system provide a comprehensive breast solution. Easy patient preparation, workflow and after-care is combined with high image quality, patient comfort and scanner efficiency. FlexTrak Mammo is available in combination with the dS Breast 16ch and/or dS Breast 7ch coil. An adjustable head rest with patient mirror and soft patient ramp are included to promote patient comfort. Height adjustment of the trolley is optional. Addition of FlexTrak Mammo to your MR system transforms it into a dedicated MR mammography solution while retaining full functionality of the Ingenia MR system.

# FlexTrak Mammo

Coverage	Bilateral
Maximum number of channels	7ch and/or 16ch
Main applications	Breast
Coil solution type	Dedicated
Data sampling	dStream digital broadband
Coil connection	dStream Interface
SENSE parallel imaging	dS SENSE enhanced parallel imaging performance



16 channel dS Breast coil



7 channel dS Breast coil



© 2015 Koninklijke Philips N.V. All rights reserved.  
Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
(Royal Philips) or their respective owners.



**How to reach us**  
Please visit [www.philips.com/healthcare](http://www.philips.com/healthcare)  
[healthcare@philips.com](mailto:healthcare@philips.com)