

Philips Avent
Niplette

Twin

## SCF152/02



# A simple solution for inverted nipples

# Effective and easy to use

Flat or inverted nipples can be a potential problem for women wishing to breastfeed. Through gentle suction the Niplette pulls the nipple out into a small plastic thimble-like cup.

### Simple solution for inverted nipples

- · Clinically proven
- Non-surgical uses gentle suction
- Simple and discreet

### When to use Nipplette

• Use with inverted or flat nipples



# **Specifications**

### **Country of origin**

England

### What is included

- Disposable breast pad: 2 pcs
- Niplette: 2 pcs

#### **Development stages**

• Stages: 0 - 6 months, Pregnancy

### Design

• Discreet design

#### Ease of use

• Easy concealing under clothing: Push syringe to release vacuum

#### Easy to clean

· Clean in hot soapy water

#### **Functions**

• Correct inverted nipples: Use gentle suction

#### Material

Breast pad: Dermatologically tested, Natural materials

# Highlights

### Non-surgical – uses gentle suction

In a matter of weeks of continued wear, the nipple will stay erect. Ideally use before becoming pregnant. However, you can also use in the first six months of pregnancy to achieve long-lasting correction. Any nipple inversion which is not congenital but has occurred recently should be immediately checked by a doctor before the Niplette is used. Use is not recommemded in the last trimester of pregnancy.

#### **Clinically proven**

The Niplette: an instrument for the non-surgical correction of inverted nipples DD McGeorge FRCS (Plast). British Journal of Plastic Surgery, December 1993.



Issue date 2016-01-23

© 2016 Koninklijke Philips N.V. All Rights reserved.

Version: 12.3.1

Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

12 NC: 0000 000 00000

www.philips.com