

# Resuscitation solutions for the operating room

Philips Internal Defibrillation Paddles

Product information



## Philips Internal Defibrillation Paddles

## We produce reusable defibrillation paddles for our Philips, Agilent and Hewlett-Packard defibrillators.

Philips reusable internal paddles are solidly constructed of top-quality materials and coverings with resilient finishes to resist cracking and protect internal circuitry from corrosion

brought on by the harshness of common cleaning methods and repeated sterilization. They are validated for steam, 100% EtO, 10/90 EtO, and Sterrad® and supported with a one-year, product replacement warranty.

Philips internal defibrillation paddles are available switched and switchless in a range of sizes from adult to infant.

As with all our devices, they are designed, tested and manufactured to the highest standards of quality and performance to enhance and support your efforts in treating patients and saving lives.

## Switchless Internal Paddles (discharge control on defibrillator)

#### Reusable, Sterilizable

Part No.	Description	Size
MI74IA	Extra Large	7.5 cm (2.95 in)
M1742A	Large	6.0 cm (2.36 in)
M1743A	Medium	4.5 cm (1.77 in)
MI744A	Small	2.8 cm (1.10 in)

## Switched Internal Paddles (discharge control on paddle grip)

#### Reusable, Sterilizable

Part No.	Description	Size
M4741A	Extra Large	7.5 cm (2.95 in)
M4742A	Large	6.0 cm (2.36 in)
M4743A	Medium	4.5 cm (1.77 in)
M4744A	Small	2.8 cm (1.10 in)

• The length of the cable from the end of the handle to the defibrillator is  $156 \pm 3$  inches for all switched and switchless paddles.

### **Internal Paddles Adapters**



# M4/40A HeartStart Internal Paddles Adapter for switchless paddles

	Adapter for	Adapter for
Defibrillator Model	<b>Switched Paddles</b>	Switchless
Paddles		
HeartStart MRx (M3535A)	Adapter Not Needed	M4740A
Heartstream/HeartStart XL (M4735A)	Adapter Not Needed	M4740A

All adapters can be cleaned with common cleaning agents — alcohol-free hand soap;
 2% gluteraldehyde solution; sodium hypochlorite (chlorine bleach), 10% in water;
 or isopropyl alcohol — for repeated use.

#### Philips—The trusted choice

- A Fortune Global 500 company, Philips is one of the world's largest medical products companies with annual revenue of over \$8 billion.
- To learn more about the HeartStart Defibrillators and Philips, call 1-800-453-6860 or visit www.philips.com/heartstart.



© 2007 Koninklijke Philips Electronics N.V. All rights are reserved.

Philips Medical Systems Nederland B.V. reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Philips Medical Systems is part of Royal Philips Electronics.

www.philips.com/heartstart medical@philips.com fax: +31 40 27 649887 Philips Medical Systems 3000 Minuteman Road Andover, MA 01810-1085

Printed in The Netherlands 4522 962 25751/861 \* SEP 2007