



Making safe connections

These neonatal blood pressure cuffs, hoses, and upgrade kits do not use Luer connectors, which the FDA and other leading organizations have said can pose a significant safety hazard. The solution was created to prevent misconnections between blood pressure devices and other medical accessories such as feeding tubes, catheters, IV tubing, and needleless IV ports. Now in both soft and vinyl materials, all Philips neonatal cuffs are DEHP-free.



NEW Soft neonatal blood pressure cuffs (with safety connectors)

Product number	Description
M1866S	Size #1 neonatal cuff
M1868S	Size #2 soft neonatal cuff
M1870S	Size #3 soft neonatal cuff
M1872S	Size #4 soft neonatal cuff
M1873S	Size #5 soft infant cuff
989803167551	Assortment - soft neonatal cuffs

Featuring soft, comfortable materials for these fragile patients, the solution is also backward compatible (with connector upgrade kits) to allow you to meet new safety requirements using older equipment.^{1,2}

**Vinyl neonatal blood pressure cuffs
(with safety connectors)**

Product number	Description
M1866B	Size #1 neonatal cuff
M1868B	Size #2 neonatal cuff
M1870B	Size #3 neonatal cuff
M1872B	Size #4 neonatal cuff
M1873B	Size #5 infant cuff
989803167541	Assortment - vinyl neonatal cuffs



Neonatal blood pressure air hoses (with safety connectors)

	Product number	Description
Air hoses	M1596C	Neonatal blood pressure air hose 1.5m
	M1597C	Neonatal blood pressure air hose 3.0m
	989803167561	Neonatal NIBP air hose
Upgrade kits	989803167511	Connector upgrade kit for M1597B
	989803167521	Connector upgrade kit for M1596B
	989803167531	Connector upgrade kit for 989803136941

1. Joint Commission on Accreditation of Healthcare Organizations (JCAHO) Sentinel Event Alert, Issue 36, April 2006; "Tubing misconnections, a persistent and potentially deadly occurrence"
2. Machinery Directive 2006/42/EC; Section 1.5.4 Errors of fitting; "Errors likely to be made when fitting or refitting certain parts which could be a source of risk must be made impossible by the design and construction of such parts"

Please visit www.philips.com/healthcarestore or call your local sales representative



© 2011 Koninklijke Philips Electronics N.V.
All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Philips Healthcare is part of Royal Philips Electronics

www.philips.com/healthcare
healthcare@philips.com

Printed in The Netherlands
4522 962 74401 * SEP 2011