



Always connected, always ready

Philips adult preconnect multifunction electrode pads

With Philips, nothing comes between you and trusted monitor/defibrillator performance. Philips preconnect electrode pads are designed to improve time to defibrillation. A single “leads-out” wire reduces cable tangling and easy-to-use packaging enables efficient deployment. These high-quality pads are validated and manufactured to precise specifications to optimize the function of your Advanced Life Support monitor/defibrillators. They are ideal for use as a preconnect pads solution by hospital, emergency medical systems, and Fire Rescue healthcare providers. These pads also enable the delivery of therapy, pacing, cardioversion, ECG monitoring and ventilation detection for Philips Q-CPR[®] measurement and feedback tool on the HeartStart MRx.

Key advantages

- Patent-pending packaging is durable and designed for fast, easy deployment
- Robust connector design delivers high-quality performance through 2,000 connection cycles
- Validated to maximize the multifunction capabilities of Philips HeartStart ALS monitor/defibrillators

PHILIPS

Philips multifunction pads means Philips quality



Durable patent-pending package design allows quick identification of pads, is easy to open, and includes graphics to guide pad placement.



Patented foam graphics, multifunction capabilities with proven gel technology, and optimal adhesive and conductive properties that best meet the needs of patients and their care providers.



Philips preconnect pads are easy to order and easy to stock, with packaging translated into multiple languages, and can be stored in a wide range of temperatures.

Robust and reliable, Philips preconnect pads are made to optimize the performance of Philips ALS monitor/defibrillators, with innovative patent-pending packaging and a connector designed to deliver high-quality performance through 2,000 connection cycles.

Philips HeartStart adult preconnect multifunction electrode pads specifications

Product number	989803166021	
Pads/case	10 pairs	
Conductor surface area	75 cm ²	
Gel surface area	102 cm ²	
Lead wire length	115 cm / 48"	
Number of hours on skin	24 hours	
Number of defibrillation shocks (EMAX)	50	
Pacing duration @ 100mA, 80 bmp	8 hours	
Temperature range for storage	0° to 43°C / 32° to 110°F	
Shelf life from date of manufacture	24 months	
Package dimensions	22.55 cm x 15.88 cm / 8.88" x 6.25"	
Connector type	Plug style	
Philips HeartStart adapter compatibility	05-10000 05-10100 05-10200 (use with M3507A)	
Philips monitor/defibrillator compatibility	M3535A	HeartStart MRx
	M3536A	HeartStart MRx
	M3536J	HeartStart MRx
	M4735A	HeartStart XL
Philips pads cable compatibility	M3508A M4763A 989803158661 M3507A (use with 05-10200)	



buy online*

Online, all the time

Our online store offers convenient ordering 24/7
www.philips.com/healthcarestore

* Available in US and other selected areas.

Please visit www.philips.com/MRxPads



© 2009 Koninklijke Philips Electronics N.V. All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Q-CPR is a registered trademark of Laerdal Medical

Philips Healthcare is part of Royal Philips Electronics

www.philips.com/healthcare
healthcare@philips.com
 fax: +31 40 27 64 887

Printed in The Netherlands
 4522 962 55871 * DEC 2009