

Value beyond measure

Philips labor and delivery supplies



Why Philips labor and delivery supplies?

Quality care

When all is said and done, the only reason to choose one brand of supplies over another is to improve the quality of care you deliver to mothers and their babies before, during, and after birth.

Philips supplies offer:

- Simplicity. Every supply we make is easy to use, giving you more time to care for mothers and their babies.
- Safety. Our new fetal spiral electrode features a protection tab that keeps the spiral tip inside during insertion.
- Accuracy. Our monitors and supplies work together to deliver accurate measurements you can trust.
- Expertise. Enjoy easy access to technical and clinical monitoring expertise — all just a phone call away.

Cost-effective care

Philips helps you manage and control the cost of labor and delivery care by making it easy to buy, use, and maintain high-quality supplies.

- Convenience. It's faster and easier to buy monitoring equipment and supplies from a single source.
- **Efficiency.** Improve the efficiency of service, maintenance, and user training.
- Support. Get all the technical, clinical, and product support you need with a single phone call.

Reliable care

Philips labor and delivery supplies deliver a total, integrated solution from a single vendor with a proven track record for accuracy and reliability in fetal and maternal monitoring.

- Compatibility. Philips supplies complement
 Philips monitors, resulting in the highest quality
 measurements.
- Accuracy. Philips world-class algorithms and high-quality supplies ensure the accuracy of our monitoring systems.
- Knowledge. The technical, product, and clinical expertise your staff needs is available from a single source.



Fetal monitoring

Non-invasive



Fetal monitoring belts


Philips broad range of soft, reusable, and disposable fetal monitoring belts includes a new Latex-free reusable model that's both water-resistant and comfortable. Available in pre-cut sizes and rolls.



Fetal monitoring paper

If you want excellent diagnostic quality and the ability to archive records long-term, you need high-quality paper. Philips fetal monitoring paper not only stands up to the rigors of labor and delivery, it maximizes the life of your Philips recorders. Our papers are mechanically and chemically compatible and provide quality traces that you depend on.

Invasive



Fetal spiral electrodes

Our accurate and reliable fetal scalp electrodes are easy to open, easy to insert, and easy to attach — for maximum fetal comfort. Compatible with all Philips, Hewlett-Packard, and Agilent fetal monitors. They require a reusable adapter cable and a disposable fetal adhesive pad electrode.

Maternal monitoring



Maternal Koala™ IUP catheters

Our disposable intrauterine sensor-tip pressure catheter is the safe, accurate way to measure IUP during contractions. Its small, soft tip is designed to avoid perforation or tissue damage. And there's no electrical risk for the mother or child.



Maternal SpO₂ sensors

Our ${\rm SpO}_2$ reusable sensors deliver instant alignment, lasting performance, accuracy, and unequalled comfort to mothers during labor and delivery. Plus, they are guaranteed for one year and last up to 18 months.



Maternal NIBP cuffs

Philips non-invasive blood pressure cuffs simplify life for clinicians, while offering moms a level of comfort that's second to none. And features like color-coded sizes, extra long cuffs, and antimicrobial coatings make this family of cuffs the right choice for labor and delivery departments.

Philips Medical Systems is part of Royal Philips Electronics

Interested?

Would you like to know more about our imaginative products? Please do not hesitate to contact us. We would be glad to hear from you.

On the web

www.medical.philips.com

Via email

medical@philips.com

By fax

+31 40 27 64 887

By mail

Philips Medical Systems Global Information Center P.O. Box 1286 5602 BG Eindhoven The Netherlands

By phone

Asia

Tel: +852 2821 5888

Europe, Middle East, Africa Tel: +49 7031 463 2254

Latin America

Tel: +55 11 2125 0764

North America

Tel: +1 800 229 6417



© 2006 Koninklijke Philips Electronics N.V. All rights are reserved.

Philips Medical Systems Nederland B.V. reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Printed in The Netherlands. 4522 962 20661/864 * NOV 2006