



PHILIPS

Prodiva 1.5T

Breast imaging solution


MR Coil

Combine your coils for **breast imaging**


A positioning mattress to position dS Flex coil that can be used in combination with the NVS Spine coil for breast imaging. The support enables bilateral breast MR imaging using dS SENSE.

Breast imaging solution


Breast imaging solution	Breast imaging mattress
Compatible coils	dS Flex Coil L&M (Combined with dS Spine coil)
Application	Breast application
Maximum number of elements	8 (2 anterior + 6 posterior elements)




Anatomy: Breast
Sequence: eThrive
Scan Time: 2.37m
Thickness: 2mm



Anatomy: Breast
Sequence: T1
Scan Time: 2.37m
Thickness: 3.5mm



Anatomy: Breast
Sequence: VISTA
Scan Time: 3.33m
Thickness: 1.8mm



Anatomy: Breast
Sequence: VISTA
Scan Time: 3.33m
Thickness: 1.8mm

© 2017 Koninklijke Philips N.V. All rights reserved.
Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

4522 991 24241 * AUG 2017



How to reach us
www.philips.com/healthcare
healthcare@philips.com

Product information
www.philips.com/mrcoils