

A woman is lying on a hospital bed, wearing a light green shirt. She has a white MR coil wrapped around her right arm. The background shows a clinical setting with a blue banner that says "Prodiva CX".

**PHILIPS**

Prodiva 1.5T

dS Flex

MR coil

## Enhance your throughput

### Flexibility to handle multiple applications

The dS Flex coil consists of two flexible elements and three sizes are available: L, M, S. The shape and size of the dS Flex coil supports a wide variety of applications and patients. The coil can be combined with the dS NVS coil to perform total neuro examinations covering head and spine.

# dS Flex coil

|                            |  |
|----------------------------|--|
| Coverage                   | 20 cm (L), 15 cm (M), 10 cm (S)  |
| Maximum number of elements | 2 elements (standalone)<br>Up to 14 elements (combined with NVS)   |
| Main applications          | dS Flex coil L: brain, hip, pelvic, posterior triangle, shoulder, pediatric (abdomen, cardiac, hip), breast<br>dS Flex coil M: shoulder, pediatric hip, hippocampus, elbow, breast, cardiac, knee<br>dS Flex coil S: TMJ, carotid artery, orbit, small joints (elbow, wrist or ankle), pediatric |
| Data sampling              | dStream digital broadband  |
| Coil connection            | dS Interface   |
| SENSE parallel imaging     | Up to 2  |
| Weight                     | ~ 0.64 kg (L), ~ 0.60 kg (M), ~ 0.58 kg (S)  |



Simplify exams with SmartSelect: Automatically detects and selects the coil and coil elements which maximize the SNR matching the area to be scanned.



This material is not for distribution in the USA.

© 2017 Koninklijke Philips N.V. All rights reserved.  
Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

4522 991 27501 \* SEP 2017



**How to reach us**  
[www.philips.com/healthcare](http://www.philips.com/healthcare)  
[healthcare@philips.com](mailto:healthcare@philips.com)

**Product information**  
[www.philips.com/mrcoils](http://www.philips.com/mrcoils)