





PHILIPS

Respironics V60 Ventilator




- Oxygen manifold kit
- Iltmanifoldsæt
- Zuurstofverdeelstuk-kit
- Hapen jakoputki
- Kit collecteur d'oxygène
- Sauerstoffverteiler-Kit
- 酸素マニホールド(オプション)
- Oksygenmanifoldsett
- Kit de tubos de oxigênio
- Kit de ramificador de oxigeno
- Syrgasförgreningsats
- Oksijen manifold kit





A



REF	Part number	O ₂ inlet connector
989805648241	1109602	DISS male
989805648741	1109881	NIST
989805650181	1113371	SIS
989805650191	1113392	DISS female (Canada)

- Your fittings may not look exactly like those shown.
- Dine beslag ligner måske ikke præcis de viste beslag.
- Uw bevestigingen zien er misschien iets anders uit dan de afbeelding.
- Omatt liittimesi saattavat näyttää hiukan erilaista kuin kuvan liittimet.
- Il se peut que vos raccords ne ressemblent pas exactement à ceux illustrés.
- Ihre Anschlüsse sehen u.U. anders aus als die abgebildeten.
- 取り付けは図示されているものと多少異なることがあります。
- Dine koblinger er ikke nødvendigvis identiske med dem som vises.
- As conexões podem não ser idênticas àquelas mostradas.
- Es posible que sus acoplamientos no tengan exactamente el mismo aspecto que los de la imagen.
- De anslutningar du har kanske inte exakt liknar de som visas.
- Parçalarınız, gösterilenlerle tamamen aynı olmayabilir.



B

Philips Healthcare is part of Royal Philips Electro

www.philips.com/healthcare
healthcare@philips.com

Respironics California, Inc.
2271 Cosmos Court
Carlsbad, CA 92011
USA
1 800 345 6443
724 387 4000

EC **REP** Respironics Deutschland GmbH
Gewerbestrasse 17
D-82211 Herrsching
Germany
+49 8 15 29 30 60

CE

© 2011 Koninklijke Philips Electronics N.V.
All rights are reserved. Reproduction or transmission in whole
or in part, in any form or by any means, electronic, mechanical
or otherwise, is prohibited without the prior written consent
of the copyright owner.

