



### Specifications

Compliance: 20mL/cmH<sub>2</sub>O  
 Connector: 15mm OD  
 Resistance: 8cmH<sub>2</sub>O/L/s  
 Breathing bag: silicone bag 1 L  
 Volume: 0-600mL  
 Weight: 230 g  
 Measurement: 300 x 100 x 43mm

Manufactured by

**GaleMed Corporation**

87, Li-gong 2<sup>nd</sup> Road

Wu-jia, I-Lan 268 Taiwan

Telephone: 886 3 990 8618

Fax: 886 3 990 8245

Website: [www.galemed.com](http://www.galemed.com)



Manufactured for

**Respironics California, Inc.**

2271 Cosmos Court

Carlsbad, CA 92011 U.S.A.

U.S.A. and Canada: 1-800-345-6443 or 724-387-4000

Europe, Africa, Middle East: +33-1-47-52-30-00

Asia Pacific: +852-3194-2280

[www.respironics.com](http://www.respironics.com)

