

Hi-Q imaging


allows you to customize your image to better identify lumen features







Hi-Q level 1



Hi-Q level 2

Choose from three high quality, optimized imaging modes with different levels of darkness and fineness of blood speckle to help differentiate between tissue and blood along the lumen border.

Available on all Refinity and Revolution catheters using Philips IVUS systems **with version 3.3 software and higher.** Hi-Q imaging modes are accessible via the "Adjust Image" menu.




Refinity ST

45 MHz rotational IVUS catheter with Hi-Q imaging

- Exceptional deliverability

 Refinity can be delivered through complex lesions with ease.
- Low profile, 5F compatibility Refinity enables radial access.
- 45 MHz high resolution imaging
 Refinity generates clear images to help facilitate accurate coronary vessel interpretation and assessment pre-, during, and post-PCI.



Ordering information

Refinity short-tip catheter	89900
Revolution catheter	89000
Revolution option (includes SpinVision device)	S5REVO 04
SpinVision sterile drapes	891BG-10 10 pack

