

**PHILIPS**

**VOLCANO**

*Verrata Plus*

Pressure guide wire



# The wire you choose, the technology you can trust



For over 20 years Philips Volcano has led the way in physiology, putting patients first through continuous innovation

## Introducing our 10th generation physiology wire: **Verrata Plus**

- Adapted from the market-leading Verrata pressure wire
- **New** in-line clip connector for improved work flow
- **New** high fidelity sensor for reliable pressure measurement
- Familiar frontline tip
- Seamlessly switch between FFR and iFR, the only resting index with over 9,000 patients studied

### Tip load (grams), n=5

\*Tip load, mean values, n=5 of each wire, data on file at Philips Volcano.



# Reliably connect and disconnect with confidence

## Robust proximal wire design

- Resists kinks
- Repels moisture

## Verrata Plus clip connector

- Multiple back up contact points
- Allows free wire rotation when connected



## Ordering information

Model	Length	Tip
10185 P	185 cm	Straight
10185J P	185 cm	J-tip
10300 P	300 cm	Straight
10300J P	300 cm	J-tip

## Technical specifications

