

dStream Head 32ch

MR coils



Advancing your neurofunctional exams

Key benefits

- Fully digital 32-element design for high SNR and fast scanning
- Provides coverage from top of the head to the brain stem and cerebellum
- Pathway for EEG leads is provided through outlet at the top of the coil

The 32-channel dStream Head coil leverages dStream technology and pushes your 3.0T system to new frontiers in neurofunctional exams. It is designed to provide high resolution volumetric imaging of the cortex, while maintaining an excellent signal-to-noise ratio (SNR) in the center of the brain. The coil supports detailed morphological and functional assessment of the brain, for instance, during pre-surgical planning. Fast scanning times are also possible through dS SENSE parallel imaging in all three measurement directions.

dStream Head 32ch coil

Systems	3.0T systems with dStream
Coverage	Anterior cranium down to the zygomatic bone
Maximum number of channels	32
Main application	Brain
Coil solution type	Dedicated
Data sampling	dStream digital broadband
Coil connection	Double FlexConnect
SENSE parallel imaging	dS SENSE enhanced parallel imaging performance
Other features	A dedicated mirror assembly for forward and rear vision is included in the package EEG procedures are supported by a dedicated pathway for the EEG leads



Axial T2w TSE Resolution: 0.4 x 0.5 x 4.0 mm Scan time: 3:21 min Ingenia 3.0T Courtesy: Hospital Q-scan, Australia



Axial 3D BrainView FLAIR Resolution: 0.6 x 0.8 x 3.0 mm Scan time: 3:46 min Ingenia 3.0T Courtesy: Hospital Q-scan, Australia



Diffusion Tensor Imaging (b800, 45 directions) Resolution: 2.0 x 2.0 x 2.2 mm Scan time: 9:59 min Achieva dStream 3.0T

© 2018 Koninklijke Philips N.V. All rights reserved. Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.



How to reach us healthcare@philips.com www.philips.com/healthcare

More information www.philips.com/mrcoils

4522 991 34451 * APR 2018