

Philips DreamStation Go
Overnight battery

Maintain therapy off the grid

Compact
Lightweight, 1.53 lbs



HH1450

CPAP therapy even 'off the grid'

For many sleep apnea patients, traveling off the power grid can create problems with maintaining therapy. Philips developed the DreamStation Go overnight battery.

Compact and lightweight

- Compact and lightweight

PHILIPS

Overnight battery

Maintain therapy off the grid Compact, Lightweight, 1.53 lbs

Specifications

Compatibility

- PAP System: DreamStation Go CPAP, AutoPAP

Service

- Warranty: 1 year

Technical information

- Weight: 696g (1.53lb)
- Output voltage: 18-24.6 VDC
- Battery technology: Lithium Ions

- Input Voltage Range: 100-240 VAC, 50/60Hz, 2.0-1.0A
- Output power: 50W
- Run time: > 8 hours
- Minimum life cycle: 500 cycles
- Recharge time: < 5 hours

Part Number

- Battery Pack: 1133281

HH1450/00

Highlights

Compact and lightweight



If you plan on traveling off the power grid, DreamStation Go's overnight battery can help you stay committed to your therapy. DreamStation Go set at 10cm of pressure displayed an average run time of 13 hours* when running off the Dreamstation Go overnight battery. Note: Battery run times may vary due to higher CPAP pressure settings, excessive mask leakage, altitude, temperature, charging devices with USB port and others. *PAP run-time test conditions: CPAP mode, pressure 10 cm H₂O, 12mm tubing, 37LPM of leak, 73.4 F room at an elevation of 1650 ft.



Issue date 2018-08-30

Version: 5.0.2

12 NC: 0000 011 33281
UPC: 6 06959 41479 2

© 2018 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com