



Philips DreamStation  
CPAP with Humidifier

Colorful, easy-to-use display

Display and track nightly use  
Designed to reduce noise  
Bluetooth enabled connectivity



HH1460/00

## Take charge with ease and confidence

With easy-to-navigate menus, a sleek, compact design and quiet operation, DreamStation makes it simple to start, customize and continue your obstructive sleep apnea (OSA) therapy.

### Performance

- Support for the long term
- Make progress every day
- Small and lightweight
- Performance and appeal


**PHILIPS**

## CPAP with Humidifier

Colorful, easy-to-use display Display and track nightly use, Designed to reduce noise, Bluetooth enabled connectivity

# Highlights

### Performance and appeal




DreamStation offers a sleek, low-profile design. Its vibrant color screen displays daily progress and features easy-to-use menus.

### Small and lightweight




DreamStation's top-mounted humidifier hose swivel and front-facing display make it easy for you to operate while lying down or sitting up in bed. And, it's small and light, making it easy to pack for travel.

### Daily Progress Feedback



Each morning you are greeted with Daily Progress Feedback - a summary of how far you've come to motivate you to stick to your therapy. The intuitive, color display shows a simple trend of their nightly hours of use, followed by a summary of the last 30 days of good nights with more than four hours of use.

### Bluetooth enabled



Built-in Bluetooth makes connectivity to the DreamMapper self-management tool easy. DreamMapper provides goal-setting tools, helpful videos and ongoing feedback to keep you motivated and engaged.

HH1460/00

# Specifications

### Product details

- Bluetooth: Built-in

### Compatibility

- For use with: DreamMapper

### Performance data

- Pressure mode: CPAP
- Flex mode: C-Flex
- Ramp time: 0-45 min (5-min increments)
- Starting ramp pressure: SmartRamp

### Operating specifications

- Operating pressure: 4 to 20 cmH<sub>2</sub>O
- Altitude compensation: Automatic



Issue date 2018-08-30

Version: 2.0.1

© 2018 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)