



**Philips InnoSpire Elegance
Compressor Nebulizer
System**

Home use table-top device

Treatment in 6-8 minutes*
*Using a 2.5ml Salbutamol dose
Nebulizer kit included



HH1300

Reliable nebulized therapy you can count on

For consistent treatment

InnoSpire Elegance provides fast, efficient and reliable aerosol drug delivery, designed for continuous treatment in the hospital and at home.

Short treatment times

- Treatment delivered in 6-8 minutes*

Stylish, robust design

- Robust design meets reliability

Proven SideStream technology

- Proven nebulizer technology

For use in hospital and home

- Designed for continuous use
- Compatible with commonly prescribed nebulizer medications

PHILIPS

Compressor Nebulizer System

Home use table-top device Treatment in 6-8 minutes*, *Using a 2.5ml Salbutamol dose, Nebulizer kit included

Highlights

For home use



Robust enough for continuous use in both a clinical and family environment, InnoSpire Elegance provides consistent and reliable nebulizer therapy.

Power and style all in one



Designed to both fit your lifestyle and cope with frequent use. You can be confident that you will receive high quality aerosol for the life of the nebulizer.

Proven nebulizer technology



As the InnoSpire Elegance pushes air through the tube, additional air is drawn in through the top of the

nebulizer to create a fast flow of aerosol medication for you to inhale. SideStream is designed to withstand frequent and repeated use and will continue to efficiently deliver your nebulizer medication until it is time for replacement.

Treatment time



InnoSpire Elegance combined with the SideStream Nebulizer, is designed to provide fast drug delivery and short treatment times.

Versatility



Capable of nebulizing a wide variety of medications that are routinely with nebulizers used to treat various respiratory conditions.

HH1300/00

Specifications

Design

- Compact: Weighs only 1.5 kg

Operating parameters

- Electrical rating/power: 120 V / 60 Hz, 1.7 Amp
- Maximum flow rate: 9.4 lpm
- Maximum pressure: 44 psi
- Mode of operation: Continuous
- Operating flow rate: 7 LPM @ 10 psi

Performance data

- MMAD: 2.9 μ m
- Fine particle fraction: 77% below 5 μ m
- Nebulization capacity: 8 ml

Product details

- Weight: 1.5 kg
- Size: 165x165x108 mm
- Noise level: 58 \pm 3 dBA
- Treatment time (6-8 minutes*): * 2.5ml Salbutamol dose

Service

- Warranty: 5 years



Issue date 2018-07-24

Version: 7.1.2

12 NC: 0000 010 99973
UPC: 3 83730 99973 9

© 2018 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

** Using a 2.5ml Salbutamol dose