

Philips Wisp
Pediatric minimal contact
nasal mask

Pack

Child size headgear
and 3 cushions

HH1038



Unique design for unique patient needs

A mask truly designed for children - from the shape of the frame to the curve of the cushion. Children can comfortably wear the mask during their normal daily activities.

Design

- Modified cushion curvature to fit an infant's face and nose
- For less intrusive mask adjustment

Performance and comfort

- Mask specifically designed for children

Child-friendly print

- Fun and child-friendly print and soft fabric materials

Open line of sight

- Read, watch TV, or wear glasses

PHILIPS


Highlights

Giraffe fabric pattern



The playful theme is to ease the child into therapy, and make him or her more comfortable with the mask.

Greater independence



Wisp's sleek design eliminates the need for a forehead pad, allowing children to comfortably wear the mask during their normal daily activities.

Innovative Leak Correction



The Leak Correction Dial lets you adjust the mask without removing it or waking the child.

Minimal contact, max comfort



The cushion, frame and headgear are designed to fit the sizing and bone structure of pediatric patients.

Soft and comfortable



All three cushion sizes fit on the same frame

Specifications

Service

- Warranty: 90-days

Design

- Compact: Performance and comfort

Specifications

- Sound pressure level: 22.1 dBA
- N° of parts: Five
- Cushion material: Silicone
- Headgear material: UBL/Polyurethane foam/Nylon/Spandex
- Tubing material: Thermoplastic polyester elastomer
- Frame material: Nylon/Spandex/Polyurethane/foam/Polyester/Acrylic
- Cleaning instructions: mask: Hand wash with mild detergent, Rinse thoroughly, air dry
- Cleaning: headgear, tubing: Hand wash with mild detergent, Rinse thoroughly, hang to dry
- Elbow/Swivel material: Amorphous Nylon
- Operating pressures: 4 -20 cmH2O
- Headgear clips material: Amorphous Nylon
- Tube wrap: Nylon/Spandex

Intent of use

- patients > 5 kg with CPAP, bi-level therapy, or non-invasive ventilation prescription



Issue date 2017-03-03

Version: 2.0.1

12 NC: 0000 011 04961
UPC: 6 06959 02673 5

© 2017 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com