



Philips SideStream  
Reusable  
Patient pack

For use with compressor system

Includes SideStream Reusable  
Includes Adult & Child Masks  
Includes tubing & mouthpiece

HH1426/00

## Treatment in 6-8 minutes\*

Proven technology makes treatment fast and easy

The SideStream Reusable Nebulizer can be used for 6 months before replacing. With its unique design, SideStream offers faster drug delivery and shorter treatment times, and can be used with most commonly prescribed respiratory medications.

### **Fast treatment times**

- Designed to draw in air

### **Consistent performance**

- Rely on consistent treatment

### **Proven technology**

- Treatment delivered in 6-8 minutes\*

**PHILIPS**

Patient pack

For use with compressor system Includes SideStream Reusable, Includes Adult & Child Masks, Includes tubing & mouthpiece

# Specifications

## Design

- Fewer parts: Easy to assemble and clean

- Venturi system: Consistent medication delivery

HH1426/00

# Highlights

## Fast, effective treatment\*




The nebulizer top, combined with an InnoSpire compressor, SideStream delivers faster aerosol drug delivery and shorter treatment times, so you can get back to your daily life.

## Effective delivery over time



Staying compliant with your breathing treatment is easier when taking your medication is consistent and reliable. SideStream nebulizers keep the driving air and the liquid drug separate for ongoing performance.

## Fast, effective treatment time



The SideStream nebulizer's venturi system boosts the airflow supplied by the compressor, delivering medication more efficiently and making your overall treatment faster.



Issue date 2018-07-27

Version: 1.0.1

© 2018 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)

\*\* Using a 2.5ml Salbutamol dose