

Philips Wisp
Headgear

Reduced size

HH1026



Click-fit headgear for easy adjustment

Elastic crown strap with integrated tube holder. Self-adjusting for a more natural fit.

Magnetic clips

- The headgear clips and mask cushion contain magnets

Design

- Personalised comfort fit

PHILIPS

Headgear
Reduced size

Specifications

Warning: Magnetic interference

- Beware with pacemakers: Read IFUs before use.

Service

- Warranty: 90 days

Design

- Comfort: Self-adjusting

Specifications

- Headgear material: UBL/Urethane foam/nylon/Lycra/nylon thermoplastic
- Cleaning: Headgear, tubing: Hand-wash with mild detergent, Rinse thoroughly, hang to dry

HH1026/00

Highlights

Magnetic clips

Contact your healthcare professional before you use this mask. Some medical devices may be affected by magnetic fields. The magnetic clips in this mask should be kept at least 50 mm (2 in.) away from any active medical device with special attention to implanted devices, such as pacemakers, defibrillators and cochlear implants. Read in-pack IFUs for more information.

Easy to adjust



Self-adjusting top strap with tube holder loop



Issue date 2018-04-24

Version: 3.1.1

12 NC: 0000 011 09307
UPC: 6 06959 40738 1

© 2018 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com