

Philips
Wisp nasal CPAP mask
with headgear

Pack
3 cushions (S/M, L, XL)
w/headgear

HH1023



Performance has a new style

Wisp's clever auto seal groove technology accommodates your every move, for a superior seal and a great night's sleep.

Magnetic clips

- The headgear clips and mask cushion contain magnets

Compact design

- Wisp comfortably fits more than 98% of users*

Open field of vision

- Read, watch TV, or wear glasses

Performance and comfort

- Provides for quick and easy fitting and adjustment.

PHILIPS

Wisp nasal CPAP mask with headgear
Pack 3 cushions (S/M, L, XL), w/headgear

Specifications

Warning: Magnetic interference

- Beware with pacemakers: Read IFU's before use.

Service

- Warranty: 90-days

Design

- Compact: Easy fitting and adjustment
- Patented design: Tip-of the nose cushion design
- Quiet exhalation: 360 degree rotating elbow
- Silicone frame: Clear unobtrusive frame
- Comfort: Elastic crown strap

Specifications

- Sound pressure level: 19 dBA
- N° of parts: Six
- Cushion material: Silicone/polycarbonate
- Headgear material: UBL/Urethane foam/nylon/Spandex/nylon thermoplastic

- Tubing material: Thermoplastic polyester elastomer
- Cleaning instructions: mask: Hand wash with mild detergent, Rinse thoroughly, air dry
- Cleaning: headgear, tubing: Hand wash with mild detergent, Rinse thoroughly, hang to dry
- Elbow/Swivel material: Polycarbonate
- Frame material: Silicone/Polycarbonate
- Headgear clips material: Acetal copolymer
- Operating pressures: 4 -20 cmH2O
- Swivel material: Silicone / Polycarbonate

Intent of use

- To be used on patients >30kg: with CPAP prescription

Compact

- Easy fitting and adjustment: Quiet exhalation

HH1023/00

Highlights

Compact design, natural fit



The 'tip-of-the-nose' cushion is designed to sit below your nose bridge, making minimal contact with your face while creating a superior seal. * Confirmation via the anthropometric database

Greater independence



Wisp's sleek design eliminates the need for a forehead pad, offering you a better field of vision.

Minimal design, superior seal



The design and headgear clips make it easier to use and clean.

Magnetic clips

Contact your healthcare professional before you use this mask. Some medical devices may be affected by magnetic fields. The magnetic clips in this mask should be kept at least 2 in. (50mm) away from any active medical device with special attention to implanted devices such as pacemakers, defibrillators and cochlear implants. Read in-pack IFU's for more information.



Issue date 2018-04-03

Version: 5.0.1

12 NC: 0000 010 94060
UPC: 6 06959 02126 6

© 2018 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

* The purchase of this medical device does not give the right to obligatory Social Security reimbursement
* Confirmation via the anthropometric database