

Philips DreamStation Go  
Reusable pollen air filter

Filters dust, pollen, & more

Replace every 6 months  
88% efficiency @ 7 - 10 micron  
Double pack

HH1453



## Reusable pollen filter

Philips reusable replacement filters help clean the air you breathe by capturing dust, pollen, mold spores, pet dander and more with an 88% efficiency at 7-10 micron.

### Reusable pollen filter

- Reusable air filter

**PHILIPS**

## Reusable pollen air filter

Filters dust, pollen, & more Replace every 6 months, 88% efficiency @ 7 - 10 micron, Double pack

# Specifications

### Compatibility

- PAP System: DreamStation Go CPAP, AutoPAP

### Service

- Warranty: 90-days

### Technical information

- Replacement: Every 6 months

- Filter material: 100% Polyester
- Filter efficiency: 88% Efficient @ 7-10 microns

### Packaging details

- Content: Includes 6 filters

### Part Number

- Reusable pollen filter: 1133743

HH1453/00

# Highlights

## Reusable air filter



Every 6 months replace your DreamStation Go reusable filter.



Issue date 2018-08-30

Version: 6.0.2

12 NC: 0000 011 33743  
UPC: 6 06959 41543 0

© 2018 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)