



Philips Amara Gel  
Headgear

Standard size



HH1220

## When it's right, it just clicks

The Amara full face mask has significantly fewer parts than leading traditional full-face masks. It requires just a single click to disassemble and reassemble the headgear.

### Minimalistic design

- 60% fewer parts

**PHILIPS**

Headgear  
Standard size

HH1220/00

# Specifications

# Highlights

## Design

- Fewer parts: Comfort

## Service

- Warranty: 90-days

## Specifications

- Headgear material: UBL/Urethane foam/nylon and Lycra
- Cleaning: headgear, tubing: Hand wash with mild detergent, Rinse thoroughly, hang to dry

## Simplicity

With only four parts, assembly and cleaning are quick and easy.



Issue date 2017-04-23

Version: 3.1.1

12 NC: 0000 010 90297

© 2017 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)