

Philips Pico
Nasal mask

Pack
w/ Small/Medium cushion

HH1032



Value-added performance

A truly appealing traditional nasal mask that makes minimal contact with the face and offers simplicity, comfort and seal.

Compact design

- Pico comfortably fits more than 98% of users*
- The lightest and smallest traditional nasal mask *

Performance and comfort

- Provides for quick and easy fitting and adjustment

PHILIPS

Nasal mask
Pack w/ Small/Medium cushion

Specifications

Light and small

- Easy fitting and adjustment: Quiet exhalation

Service

- Warranty: 90-days

Specifications

- Operating pressure: 4 -20 cm H₂O
- Sound pressure level: 19 dBA
- N° of parts: Five
- Cushion material: Silicone/polycarbonate
- Headgear material: UBL/Urethane foam/nylon/

- Lycra
- Cleaning instructions: mask: Hand wash with mild detergent, Rinse thoroughly, air dry
- Cleaning: headgear, tubing: Hand wash with mild detergent, Rinse thoroughly, hang to dry
- Frame/Elbow/Swivel material: Polycarbonate
- Talon clips material: Thermoplastic polyester elastomer

Intent of use

- To be used on patients >30kg: with CPAP prescription

HH1032/01

Highlights

Compact design for a natural f




The 'tip-of-the-nose' cushion is designed to sit below your nose bridge, making minimal contact with your face while creating a superior seal. * Data on file, Internal Marketing Survey 2013

Minimal design, superior seal



The design and headgear clips make it easier to use and clean.

The lightest, smallest choice



Pico makes minimal contact with the face, providing the freedom to watch television or read while wearing it. * Data on file, Internal Marketing Survey 2013



Issue date 2017-03-06

Version: 1.0.1

12 NC: 0000 011 04921
UPC: 6 06959 00893 9

© 2017 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com