



Philips DreamWear
Frame

Replacement Part
Small frame



HH1116

Like wearing nothing at all.*

With the freedom and comfort of DreamWear, you can sleep how you want.

Innovative design

- Air tube connection at the top of the head
- Soft, flexible frame

Easy to use

- Features just for you

PHILIPS

Specifications

Innovative design

- Comfort and performance: Soft fabric straps

Product details

- Materials: Polycarbonate, Silicone
- Warranty: 90 days

Cleaning: Frame


- Hand-wash with mild detergent: Rinse thoroughly, air-dry¹

Replacement

- Inspect daily for wear: Replace at first sign of wear

Highlights

Easy-to-use frame



Easy to use, choose one of the three frame sizes that fits you best. All of the frame sizes fit all of the cushion sizes.

Innovative frame design



The soft, flexible frame is comfortable to wear at night. Users have said that it feels like wearing nothing at all.*

Freedom of movement



By moving the tube location to the top of the head, DreamWear offers more freedom of movement during the night. Users said that they didn't need to choose a sleeping position when wearing DreamWear.*



Issue date 2018-04-03

Version: 10.0.1

12 NC: 0000 011 16745
UPC: 6 06959 40204 1

© 2018 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

* *2017 Philips User Preference Questionnaire

** Data analysis of Aug 2017 Patient Preference trial where n = 85
and prescribed masks include Resmed Airfit F10 and F20,
Respironics Amara View and F&P Simplus; based on 10 days of use.

*¹ Please refer to the Instructions For Use for more detailed cleaning
instructions.