



Philips
Amara View full-face mask

Fit pack — 3 cushion sizes
1 modular frame

Enjoy the view

With its innovative under-the-nose design, Amara View offers an open field of vision and prevents red marks, discomfort or irritation on the bridge of the nose.



HH1231

Open field of vision

- Widest field of vision of all leading full-face masks

Innovative design

- Headgear with crown strap
- Four-part modular design
- Innovative design

The smallest and lightest choice

- The smallest* and lightest** of all leading full-face masks

PHILIPS

Amara View full-face mask
Fit pack — 3 cushion sizes 1 modular frame

Highlights

Smallest and lightest choice



Smaller and lighter than traditional full-face masks, Amara View covers less of the patient's face than comparable masks and eliminates the bulky cushion and frame in front of the eyes.

Unobstructed view



It's easy to wear glasses, read, watch TV and use a computer or tablet before falling asleep.

Headgear with crown strap



Headgear with soft and comfortable fabric straps. The crown strap provides quick and easy fitting adjustment.

Four-part modular design



Modular design for quick assembly and cleaning, with interchangeable cushions available in four sizes.

Under-the-nose cushion



Prevents red marks, discomfort or irritation on the bridge of the nose

HH1231/00

Specifications

Service

- Warranty: 90 days

Innovative design

- Wide field of vision: Wear glasses, read or watch TV

Contraindications:

- Impaired cough reflex: impaired cardiac sphincter
- Recent eye surgery, dry eyes; hiatal hernia, major reflux

Intent of use

- To be used on patients >30 kg (>66 lbs): with CPAP prescription

Specifications

- Operating pressure: 5–30 cmH₂O
- Sound pressure level: 24.2 dBA
- No. of parts: Four
- Cushion material: Silicone
- Headgear material: UBL/mesh blend
- Cleaning instructions: mask: Rinse thoroughly, air dry
- Cleaning: headgear, tubing: Hand-wash with mild detergent, Rinse thoroughly, hang to dry
- Headgear clips/Tubing material: Thermoplastic polyester elastomer
- Quick-release feature material: Polyoxymethylene (Acetal)
- Frame/Elbow/Swivel material: Polycarbonate



Issue date 2018-03-14

Version: 4.0.1

12 NC: 0000 010 90670
UPC: 6 06959 00928 8

© 2018 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

* Internal testing: data on file

** Excludes Amara View's quick-release tube

* The purchase of this medical device does not give the right to obligatory Social Security reimbursement