



Philips DreamWear
Minimal contact nasal mask

Pack

and 4 cushion sizes
1 frame size Medium



HH1111

Like wearing nothing at all*

With the freedom and comfort of DreamWear, you can sleep how you want *2015
Philips User Preference Questionnaire

The best of pillows and nasal masks

- Comfort and performance

Comfortable design

- Slim-line headgear
- Soft, flexible frame

Easy to use

- Comfort and soft

Innovative design

- Under-the-nose cushion
- Air tube connection at the top of the head
- Four-part modular design

Open field of vision

- Open design offers an extremely wide field of vision

PHILIPS

Minimal contact nasal mask
Pack and 4 cushion sizes, 1 frame size Medium

HH1111/00


Highlights

Innovative under the nose mask



Prevents red marks, discomfort or irritation in the nostrils or on the nose bridge*

Freedom of movement



By moving the tube location to the top of the head, DreamWear offers more freedom of movement during the night. Users stated they didn't need to choose a sleep position when wearing DreamWear.*

Four-part modular design



Modular design for quick assembly and cleaning, with interchangeable cushions available in four sizes.

Slim-line headgear



DreamWear combines comfort with aesthetics: Not only will you barely feel the mask, it is also subtle and unobtrusive.

Soft fabric wraps



Provides a soft touch on your cheek

Soft, flexible frame



The flexible material through which the air flows is applied gently to the skin and is pleasant to wear.

Under-the-nose cushion



Prevents red marks and irritation on the bridge of the nose that many patients experience while wearing a traditional full-face mask

Unobstructed view



It's easy to wear glasses, read, watch TV, and use a computer or tablet before falling asleep.

Specifications

Innovative design

- Open field of vision: Wear glasses, read or watch TV

Service

- Warranty: 90-days

Comfort

- Soft fabric wraps: Soft interface to the cheek
- Soft, flexible frame: Gives a soft touch on the face
- Under-the-nose cushion: Prevents red marks in nostrils

Intent of use

- To be used on patients >30kg: with CPAP prescription

Specifications

- Cleaning instructions: mask: Hand wash with mild detergent, Rinse thoroughly, air dry
- Cleaning: headgear, tubing: Hand wash with mild detergent, Rinse thoroughly, hang to dry
- Cushion material: Silicone
- Elbow/Swivel material: Polycarbonate
- Frame material: Polycarbonate
- Headgear material: Polyurethane foam, nylon/Spandex
- N° of parts: Five
- Operating pressure: 4 -20 cmH₂O
- Sound pressure level: 18 dBA
- Fabric sleeves material: Polyester/Spandex



Issue date 2018-04-03

Version: 3.0.1

12 NC: 0000 011 16720
UPC: 6 06959 02340 6

© 2018 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com