



Philips DreamStation
Travel kit for
DreamStation CPAP
device

Case color: light gray
Including 65W Power Supply
with USB port



HH1423

Ideal for travelling

Compact travel bag for DreamStation

Gives patients a convenient, compact, and discreet way to travel with DreamStation

Ideal combination of protection and portability

- Benefit from your CPAP therapy on the move

Ideal for travel

- Dedicated compartments

Including 65W Power Supply with USB Port

- with USB Port

PHILIPS

Travel kit for DreamStation CPAP device
Case color: light gray Including 65W Power Supply, with USB port

HH1423/00

Specifications

Service

- Warranty: 90-days

Travel

- Dedicated to DreamStation: Travel case with 65W power supply with USB port holds mask, headgear, tubing and device

Highlights

Lightweight power supply



The case includes a 65W power supply which is smaller than the power supply used for the DreamStation machine and humidifier. The power supply has a built-in USB port for charging cell phones and other electronics while using the DreamStation machine.

Compact and easy to use



A molded insert liner helps the patient thoroughly pack the PAP system components: mask and tubing, 65W power supply, DreamStation WITHOUT humidifier.

Easy to care and sturdy



Lightweight, hard-fabric case. Removable shoulder strap + trolley back strap to easily slide onto luggage handles. Secure zippered closures help facilitate airport screening; and the travel case easily fits in airplane overhead bins and under airline seats.



Issue date 2017-10-04

Version: 5.0.2

12 NC: 000E U11 20135
UPC: 6 06959 02405 2

© 2017 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

* The purchase of this medical device does not give the right to obligatory Social Security reimbursement

* These comfort accessories are optional. They can be purchased in addition, but are not part of the benefit covered by the Health Insurance Scheme.