

Philips Wisp Frame

# Soft and pliable fabric frame

Reversible fabric frame in soft suede or silky sateen. For a more natural feeling

### with fabric frame

## Fabric frameChoice of soft suedette or silky sateen



HH1000



### Specifications

#### **Service**

• Warranty: 90 days

#### Design

- Compact: Easy fitting and adjustment
  Fabric frame: Soft suede or silky sateen
- Comfort: Elastic crown strap
- Quiet exhalation: 360-degree rotating elbow
- · Patented design: Tip-of the nose cushion design

#### **Specifications**

• Frame material: Nylon/Lycra/Urethane foam

#### HH1000/00

### Highlights

#### Soft and comfortable



You told us that you wanted a better choice of mask, something that combined the performance and comfort of a nasal mask with the aesthetic elegance of a pillows mask.



Issue date 2018-04-03 © 2018 Koninklijke Philips N.V. All Rights reserved. Version: 5.3.5 Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

12 NC: 0000 010 94081 UPC: 6 06959 00688 1

www.philips.com